RAPIDVUE[®]

● *hCG test*

CLIA Complexity: WAIVED

INTENDED USE

The RapidVue hCG test is a sensitive immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine for the early detection of pregnancy.

SUMMARY AND EXPLANATION

Human chorionic gonadotropin is a hormone produced by the placenta shortly after implantation. Since hCG is present in the urine of pregnant women, it is an excellent marker for confirming pregnancy.

The RapidVue hCG test uses a monoclonal antibody specific to the beta subunit of hCG in a novel single-step technology to accurately detect hCG as early as 2 or 3 days before the expected menses.

PRINCIPLE OF THE TEST

The Test Strip is briefly dipped in urine. The sample migrates by capillary action through a pad containing a monoclonal anti-hCG antibody; if the sample contains hCG, a complex is formed with the anti-hCG antibody. The complex will continue to migrate and be bound by a polyclonal anti-hCG antibody spotted on the test line.

If hCG is present in the sample at a level of 25 mIU/mL or greater, a pink Test (T) Line along with a blue Control (C) Line will appear. If hCG is not present in the sample, or present at low levels, only a blue Control Line will appear.

REAGENTS AND MATERIALS SUPPLIED

- 25 Test Strips
 - ▶ Each strip contains murine monoclonal antibody and caprine polyclonal antibody to hCG
- 1 Package Insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer or watch that measures minutes and seconds
- Sample collection containers

MATERIALS RECOMMENDED BUT NOT PROVIDED

External hCG controls traceable to an appropriate WHO Standard, preferably WHO 3rd IS 75/537 hCG Standard (see Quality Control section).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- **Use Test Strip within 15 minutes of removal from pouch.**
- Hold the Test Strip only by the top (plastic-covered end).
- Do not bend Test Strips.
- Do not dip the Test Strip past the line marked on the Strip indicated by 3 arrows.
- **Do not use kit contents after the expiration date printed on the outside of the kit.**
- Dispose of containers and used contents in accordance with Federal, State and Local requirements.³
- To obtain accurate results, you must follow the Package Insert instructions.

KIT STORAGE AND STABILITY

Store kit at room temperature 59–86°F (15–30°C) out of direct sunlight. Kit contents are stable until the expiration date printed on the outer box.

SPECIMEN COLLECTION AND STORAGE

Collect urine samples in a clean container. Urine samples may be kept at room temperature for 8 hours or stored at 2–8°C for up to 3 days. Do not freeze samples. Do not mix stored samples.

QUALITY CONTROL

Built-in Control Features

The RapidVue hCG test contains built-in control features. The development of the blue Control Line indicates that proper volume of fluid was absorbed into the Test Strip and capillary flow occurred. If the Control Line does not develop within 3 minutes, the test result is invalid. If background color appears which interferes with your ability to read the test result, your result may be invalid. In this case, review the procedure and repeat test with a new Test Strip.

Positive and Negative Quality Control

Good laboratory practice recommends the use of external controls to assure that the assay is performing properly. The external controls used by your laboratory must be traceable to an appropriate WHO Standard. For this purpose we recommend the Quidel hCG Controls (Quidel Catalog Number 00272). It is recommended that the Quidel hCG Controls be tested once for each 25 tests and as otherwise required by your laboratory's standard quality control procedures.


If the external controls do not perform as expected, do not use the test results. Repeat the test or contact Quidel Technical Support.

TEST PROCEDURE


NOTE: Several tests may be run at one time.

Set up a sufficient number of appropriately labeled Test Strips.

- Open at notch. Remove a RapidVue hCG test Strip from the foil pouch.


- Dip the Test Strip only up to the line indicated on the strip by 3 arrows in the urine sample for 10 SECONDS.

- REMOVE the Test Strip from the urine sample and place it on a clean, flat, dry surface. The Test Strip should not be jarred or manipulated until the assay is ready for interpretation


- Read results at THREE MINUTES. Some positive results may appear earlier.

INTERPRETATION OF RESULTS


POSITIVE:

The appearance of any pink-to-purple Test Line along with a blue Control Line.


NEGATIVE:

The appearance of the blue Control Line only.

NO RESULT:

A blue Control Line should always appear. If no blue Control Line appears, the test is invalid, and the sample must be re-tested.

If a negative result is obtained but pregnancy is suspected, another sample should be collected after 48–72 hours and tested.

LIMITATIONS

The contents of this kit are for use in the qualitative detection of hCG in urine only.

Conditions other than normal pregnancy may be associated with detectable hCG, including, for example, ectopic pregnancy or molar pregnancy.^{1,2} hCG may remain detectable for a few days to several weeks after delivery, abortion, or hCG injections.⁴ Natural termination occurs in 31% of pregnancies overall and 22% of clinically unrecognized pregnancies.⁵ Abnormal pregnancies cannot be diagnosed by qualitative hCG results. The above conditions should be ruled out when diagnosing pregnancy.

NOTE: A SAMPLE WITH A LOW LEVEL OF hCG MAY SHOW COLOR DEVELOPMENT OVER TIME.

Test results must always be evaluated with other data available to the physician.

EXPECTED VALUES

Samples containing as low as 25 mIU/mL (calibrated against the WHO 3rd IS 75/537) hCG will yield positive results when tested with the RapidVue hCG test.⁶ In normal pregnancy, hCG can be detected as early as 6 days following conception with concentrations doubling every 32 to 48 hours. For some patients, an hCG level of 25 mIU/mL can be detected as early as two to three days before expected menses. In general, an hCG level of 50 mIU/mL is found approximately 10 days after conception.⁷ By the first day of missed menses, concentrations often exceed 100 mIU/mL.

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity, Specificity, and Accuracy

The RapidVue hCG test will yield positive results with samples containing hCG levels as low as 25 mIU/mL.

The RapidVue hCG test was used to assay urine samples from patients presenting for pregnancy testing. Positive and negative results were compared to results obtained with two different commercially available, visually interpreted tests for hCG. A commercially available quantitative radioimmunoassay (RIA) was used to resolve any discrepant results.

Of the 740 urine samples evaluated in a multi-center field study, 232 samples tested positive and 506 samples tested negative by both the RapidVue hCG test and the comparative test methods. An agreement of 99.7% was determined for these samples.

		Correlation Comparative Test	
		+	-
RapidVue hCG	+	232	1
	-	1	506

Sensitivity: $232/233 = 99.6\%$

Specificity: $506/507 = 99.8\%$

Accuracy: $738/740 = 99.7\%$

Physician's Office Laboratory (POL) Studies

An evaluation of the RapidVue hCG test was conducted at three Physician's Offices using a panel of coded specimens. Testing was performed by physician's office personnel with diverse educational backgrounds and work experience at different locations. The proficiency panel contained negative, low positive, moderate positive and high positive specimens. Each level was tested in replicates of five at each site over a period of three days.

Not to be used for performing assay. Refer to most current package insert accompanying your test kit.

The results at each site agreed with the expected results. No significant differences were observed within run (5 replicates), between runs (three different assay days), or between sites (three POL sites).

Cross-Reactivity

hTSH, hLH, and hFSH were tested in the RapidVue hCG test at levels ranging from 1000 µIU/mL to 1000 mIU/mL and did not affect the expected results.

Interference Testing

The following chemical and biological compounds were tested in the RapidVue hCG test and did not affect the expected results.

Table 1

Chemical Analytes	Concentration
Acetaminophen	20 mg/dL
Ascorbic Acid	20 mg/dL
Caffeine	20 mg/dL
Gentisic Acid	20 mg/dL
Salicylic Acid	20 mg/dL
EDTA	80 mg/dL
Cannabinol	10 mg/dL
Methadone	10 mg/dL
Methanol	10.0%
Acetoacetic Acid	2000 mg/dL
β-Hydroxybutyrate	2000 mg/dL
Ephedrine	20 mg/dL
Phenylpropanolamine	20 mg/dL
Phenothiazine	20 mg/dL
Acetylsalicylic Acid	20 mg/dL
Benzoyllecgonine (cocaine metabolite)	10 mg/dL
Ethanol	1.0%
DMSO	5.0%
Uric Acid	20 mg/dL
Urine Analytes	Concentration
Albumin (serum)	2000 mg/dL
Bilirubin	1000 µg/dL
Hemoglobin	1000 µg/dL
Glucose	2000 mg/dL
Urine pH	5–9

Hormones	Concentration
hLH	1000 mIU/mL
hFSH	1000 mIU/mL
hTSH	1000 μ IU/mL
Estriol 17-beta	1400 μ g/mL
Pregnanediol	1500 μ g/mL

Bacteria	Concentration
E. coli	10 ⁸ CFU/mL
Group B <i>Streptococcus</i>	2.5 x 10 ⁸ CFU/mL
<i>Chlamydia trachomatis</i>	10 ⁴ IFU/mL

ASSISTANCE

If you have any questions regarding the use of this product, please call Quidel's Technical Support number, (800) 874-1517 (toll-free in the U.S.A.) or (858) 552-1100, Monday through Friday, between 7:00 a.m. and 5:00 p.m. Pacific Time, U.S.A. If outside the United States, contact your local distributor or technicalsupport@quidel.com.

REFERENCES

1. Saxena BB. In *Endocrinology of Pregnancy*, 3rd ed., F Fuchs, A Klopper, Eds., Harper and Row, Philadelphia, PA, 1983, pp 50-72.
2. Krieg AF, In *Clinical Diagnosis and Management by Laboratory Methods*, Vol. 1, 16th ed., JB Henry, Ed., W.B. Saunders Co., Philadelphia, 1979, pp 680-692.
3. *Recommendations for the Prevention of HIV Transmission in Health Care Settings*, Morbidity and Mortality Weekly Report, Centers for Disease Control, August 21, 1987.
4. Steier JA, Bergsjö P, Myking OL. *Obstet Gynecol* 64: 391-394 (1984).
5. Wilcox AJ, Weinberg CR, O'Connor JF, Baird DD, Schlatterer JP, Canfield RE, Armstrong EG, Nisula BC. *Incidence of Early Loss of Pregnancy*, *N Eng J Med* 319: 189-194 (1988).
6. Twenty-Sixth Report of WHO Expert Committee on Biological Standardization. *WHO Tech Report Series*, No. 565, 1975.
7. McCready J, Braunstein GD, Helm D, Wade ME, *Clin Chem* 24: 1958-1961, (1978).

RapidVue hCG test covered by U.S. Patent Nos. 4,943,522, 5,766,961 and 5,770,460; European Patent Nos. 0 296 724 and 0 566 695; other patents pending.

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Quidel Corporation
Worldwide Headquarters
10165 McKellar Court
San Diego, CA 92121 USA
www.quidel.com

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LITERATUR

1. Saxena BB. In *Endocrinology of Pregnancy*, 3rd ed., F Fuchs, A Klopper, Eds., Harper and Row, Philadelphia, PA, 1983, pp 50-72.
2. Krieg AF, In *Clinical Diagnosis and Management by Laboratory Methods*, Vol. 1, 16th ed., JB Henry, Ed., W.B. Saunders Co., Philadelphia, 1979, pp 680-692.
3. *Recommendations for the Prevention of HIV Transmission in Health Care Settings*, Morbidity and Mortality Weekly Report, Centers for Disease Control, August 21, 1987.
4. Steier JA, Bergsjö P, Myking OL. *Obstet Gynecol* 64: 391-394 (1984).
5. Wilcox AJ, Weinberg CR, O'Connor JF, Baird DD, Schlatterer JP, Canfield RE, Armstrong EG, Nisula BC. *Incidence of Early Loss of Pregnancy*, N Eng J Med 319: 189-194 (1988).
6. Twenty-Sixth Report of WHO Expert Committee on Biological Standardization. *WHO Tech Report Series*, No. 565, 1975.
7. McCready J, Braunstein GD, Helm D, Wade ME, *Clin Chem* 24: 1958-1961, (1978).

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REFERENZE

1. Saxena BB. In *Endocrinology of Pregnancy*, 3rd ed., F Fuchs, A Klopper, Eds., Harper and Row, Philadelphia, PA, 1983, pp 50-72.
2. Krieg AF, In *Clinical Diagnosis and Management by Laboratory Methods*, Vol. 1, 16th ed., JB Henry, Ed., W.B. Saunders Co., Philadelphia, 1979, pp 680-692.
3. *Recommendations for the Prevention of HIV Transmission in Health Care Settings*, Morbidity and Mortality Weekly Report, Centers for Disease Control, August 21, 1987.
4. Steier JA, Bergsjø P, Myking OL. *Obstet Gynecol* 64: 391-394 (1984).
5. Wilcox AJ, Weinberg CR, O'Connor JF, Baird DD, Schlatterer JP, Canfield RE, Armstrong EG, Nisula BC. *Incidence of Early Loss of Pregnancy*, N Eng J Med 319: 189-194 (1988).
6. Twenty-Sixth Report of WHO Expert Committee on Biological Standardization. *WHO Tech Report Series*, No. 565, 1975.
7. McCready J, Braunstein GD, Helm D, Wade ME, *Clin Chem* 24: 1958-1961, (1978).

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