

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17

*For use under the Emergency Use Authorization  
(EUA) only  
For in vitro diagnostic use  
Rx Only*

**Lyra® SARS-CoV-2 Assay  
Instructions for Use**

*For the qualitative detection of human coronavirus SARS-CoV-2 viral RNA extracted from nasal, nasopharyngeal and oropharyngeal swab specimens.*

18	<b>Contents</b>	
19	Intended Use .....	4
20	Summary and Explanation.....	4
21	Principle of the Procedure.....	4
22	Materials Provided .....	5
23	Materials Required But Not Provided .....	6
24	Warnings and Precautions.....	6
25	Storage and Handling of Kit Reagents .....	7
26	Indications of Instability or Deterioration of Reagents .....	7
27	Specimen Collection, Storage and Handling .....	7
28	Nucleic Acid Extracts Storage .....	8
29	bioMérieux NucliSENS easyMAG Nucleic Acid Extraction Programming Instructions.....	8
30	Assay Procedure .....	10
31	Master Mix Rehydration Procedure .....	10
32	RT-PCR Set-up Procedure .....	10
33	Quality Control .....	10
34	Interpretation of Results from Patient Specimens.....	11
35	CLINICAL PERFORMANCE .....	12
36	ANALYTICAL PERFORMANCE .....	13
37	Level of Detection.....	13
38	Analytical Reactivity (Inclusivity) .....	17
39	Analytical Specificity (Cross-Reactivity).....	17
40	Limitations .....	19
41	Conditions of Authorization for the Labs .....	20
42	Customer and Technical Assistance .....	21
43	Intellectual Property.....	21

44	References .....	21
45	APPENDIX .....	22
46	Applied Biosystems 7500 Fast Dx Programming Instructions .....	22
47	Applied Biosystems 7500 Fast Dx Thermocycler Test Procedure .....	23
48	Applied Biosystems 7500 Standard Programming Instructions .....	24
49	Applied Biosystems 7500 Standard Thermocycler Test Procedure .....	26
50	Bio-Rad CFX96 Touch Thermocycler Programming Procedure .....	27
51	Bio-Rad CFX96 Touch Thermocycler Test Procedure .....	29
52	Qiagen Rotor-Gene Q Programming Instructions .....	29
53	Qiagen Rotor-Gene Q Test Run .....	31
54	Roche LightCycler 480 Instrument II Programming Instructions .....	31
55	Creating a LightCycler 480 II Assay Run Template .....	31
56	Creating a Roche LightCycler 480 II Assay Test Procedure.....	33
57	Thermo Fisher QuantStudio 7 Pro Programming Instructions .....	33
58	Thermo Fisher QuantStudio 7 Pro Test Run Programming Instructions:.....	33
59	Creating a Thermo Fisher QuantStudio 7 Pro Test Procedure .....	35
60	GLOSSARY .....	38
61		
62		

63 **Intended Use**

64 The Lyra SARS-CoV-2 Assay is a real-time RT-PCR assay intended for the qualitative detection of nucleic  
 65 acid from SARS-CoV-2 in nasal, nasopharyngeal (NP), or oropharyngeal (OP) swab specimens from  
 66 patients suspected of COVID-19 by their healthcare provider. Testing is limited to laboratories certified  
 67 under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, to perform high  
 68 complexity tests.

69  
 70 Results are for the identification of SARS-CoV-2 RNA. The SARS-CoV-2 is generally detectable in upper  
 71 respiratory specimens during the acute phase of infection. Positive results are indicative of the presence  
 72 of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary  
 73 to determine patient infection status. Positive results do not rule out bacterial infection or co-infection  
 74 with other viruses. Laboratories within the United States and its territories are required to report all  
 75 positive results to the appropriate public health authorities.

76  
 77 Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient  
 78 management decisions. Negative results must be combined with clinical observations, patient history, and  
 79 epidemiological information.

80  
 81 The Lyra SARS-CoV-2 Assay is intended for use by qualified and trained clinical laboratory personnel  
 82 specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures.  
 83 The Lyra SARS-CoV-2 Assay is only for use under the Food and Drug Administration’s Emergency Use  
 84 Authorization.

85  
 86 **Summary and Explanation**

87 SARS-CoV-2, also known as the COVID-19 virus, was first identified in Wuhan, Hubei Province, China  
 88 December 2019. This virus, as with the novel coronavirus SARS-1 and MERS, is thought to have originated  
 89 in bats, however the SARS-CoV-2 may have had an intermediary host such as pangolins, pigs or civets.<sup>1</sup> By  
 90 the start of March 2020, human infection has spread to over 74 countries, infected over 92,000 people and  
 91 has killed over 3100 people.<sup>1</sup> On March 11, the WHO had declared the SARS-CoV-2 as a global pandemic.

92  
 93 The median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12  
 94 days of infection.<sup>2</sup> The symptoms of COVID-19 are similar to other viral respiratory diseases and include  
 95 fever, cough and shortness of breath.<sup>3</sup>

96  
 97 The Lyra SARS-CoV-2 Assay has been designed to specifically detect SARS-CoV-2 RNA.

98  
 99 **Principle of the Procedure**

100 The Lyra SARS-CoV-2 Assay detects SARS-CoV-2 viral RNA that has been extracted from a patient sample  
 101 using either the bioMérieux NucliSENS® easyMAG® system or EMAG® system. A multiplex real-time RT-PCR  
 102 reaction is carried out under optimized conditions in a single tube generating amplicons for the targeted  
 103 virus (if present) and the Process Control (PRC) present in the sample. This reaction is performed utilizing  
 104 one of six thermocyclers: Applied Biosystems® 7500 Fast Dx, Applied Biosystems 7500 Standard, Roche  
 105 LightCycler® 480, Qiagen® Rotor-Gene Q, Bio-Rad® CFX96 Touch™, Thermo Fisher QuantStudio™ 7 Pro.  
 106 Identification of the SARS-CoV-2 virus occurs by the use of target specific primers and fluorescent-labeled  
 107 probes that hybridize to a conserved region of the non-structural polyprotein of the SARS-CoV-2 virus.

108

<b>Table 1. Lyra® SARS-CoV-2 Assay Probe Labels</b>	
<b>Target</b>	<b>Dye</b>
Non-structural polyprotein (pp1ab)	FAM
Process Control (PRC)	Quasar® 670

110 The following is a summary of the procedure:

111

112 1. **Sample Collection:** Obtain nasopharyngeal, oropharyngeal, or nasal swabs using standard techniques from  
113 symptomatic patients. These specimens are transported, stored, and processed according to established  
114 laboratory procedures.

115

116 2. **Nucleic Acid Extraction:** Extract nucleic acids from the specimens with the NucliSENS easyMAG or EMAG  
117 Systems following the manufacturer’s instructions and using the appropriate reagents (See **Materials**  
118 **Required but Not Provided**).

119

120 Prior to the extraction procedure add 20 µL of the Process Control (PRC) to each 180 µL aliquot of specimen  
121 or controls. The PRC serves to monitor inhibitors in the extracted specimen, assures that adequate  
122 amplification has taken place and confirms that the nucleic acid extraction was sufficient.

123

124 3. **Rehydration of Master Mix:** Rehydrate the lyophilized Master Mix using 135µL of Rehydration Solution.  
125 The Master Mix contains oligonucleotide primers, fluorophore and quencher-labeled probes targeting  
126 conserved regions of the SARS-CoV-2 as well as the process control sequence. The probes are dual labeled  
127 with a reporter dye attached to the 5’ end and a quencher attached to the 3’ end. The rehydrated Master  
128 Mix is sufficient for eight reactions.

129

130 4. **Nucleic Acid Amplification and Detection:** Add 15 µL of the rehydrated Master Mix to each plate well  
131 (Applied Biosystems 7500 Fast Dx, Applied Biosystems 7500 Standard, the Roche LightCycler 480) or tube  
132 (Qiagen Rotor-Gene Q). 5 µL of extracted nucleic acids (specimen with PRC) is then added to the plate well  
133 or tube. Place the plate or tube into the appropriate instrument.

134

135 Once the reaction plate or tubes are added to the instrument, the assay protocol is initiated. This protocol  
136 initiates reverse transcription of the RNA targets generating complementary DNA, and the subsequent  
137 amplification of the target sequences occurs. The Lyra SARS-CoV-2 Assay is based on TaqMan® chemistry,  
138 and uses an enzyme with reverse transcriptase, DNA polymerase, and 5’-3’ exonuclease activities. During  
139 DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating  
140 the quencher dye from the reporter dye. This step generates an increase in fluorescent signal upon  
141 excitation by a light source of the appropriate wavelength. With each cycle, additional dye molecules are  
142 separated from their quenchers resulting in additional signal. If sufficient fluorescence is achieved, the  
143 sample is reported as positive for the detected target sequence.

144

145 **Materials Provided**

146

SKU # CE-M120

147

#	Component	Quantity
①	<b>Rehydration Solution</b> Part M5003	1 vial/kit 1.9 mL
②	<b>Lyra SARS-CoV-2 Master Mix</b> Part M5150 Lyophilized Contents: DNA polymerase enzyme with reverse transcriptase activity Oligonucleotide primer pairs; Oligonucleotide probes dNTPs (dATP, dCTP, dGTP, dUTP, dTTP) Stabilizers	12 vials/kit, 8 reactions/vial
<b>CONTROL</b>	<b>Process Control</b> Part M5273	1 vial/kit 2.0 mL
<b>CONTROL+</b>	<b>Positive Control</b> containing SARS-CoV-2 Synthetic RNA, Part M5153	1 vial/kit 1.0 mL

<b>Table 2. Detection Kit (96 Reactions) – Store at 2°C to 8°C</b>		
<b>#</b>	<b>Component</b>	<b>Quantity</b>
<b>CONTROL-</b>	<b>Negative Control</b> Part M5275	1 vial/kit 2.0 mL

148

- 149 • Lyra SARS-CoV-2 Assay Instructions for Use

150

### 151 **Materials Required But Not Provided**

- 152 • Micropipettors (range between 1 to 10 µL and 100 to 1000 µL)

- 153 • Non-aerosol pipette tips

- 154 • Applied Biosystems 7500 Fast Dx, software version 1.4

- 155 • Applied Biosystems 7500 Standard, software version 2.0.6

- 156 • Roche LightCycler 480 Instrument II, software version 1.5.0.39

- 157 • Qiagen Rotor-Gene Q, software version 2.0.2.4

- 158 • Bio-Rad CFX96 Touch, software version 3.1

- 159 • Thermo Fisher QuantStudio 7 Pro, software version 2.0

- 160 • 96 well PCR plate #:

- 161     – Applied Biosystems 7500 Fast Dx: 4344906

- 162     – Applied Biosystems 7500 Standard: N8010560

- 163     – Roche LightCycler 480: 04729692001, foil included

- 164     – Bio-Rad CFX96 Touch: HSP9631, seals MSB1001

- 165     – Thermo Fisher QuantStudio 7 Pro: 4483354

- 166 • Optical plate films

- 167 • Qiagen Rotor-Disc

- 168 • Qiagen Rotor-Disc Heat Sealing Film

- 169 • Plate centrifuge for 96 well plate

- 170 • bioMérieux NucliSENS easyMAG software version 2.0

- 171 • bioMérieux EMAG software version 2.0

- 172 • bioMérieux NucliSENS easyMAG Buffers 1, 2, 3

- 173 • bioMérieux NucliSENS easyMAG Lysis Buffer

- 174 • bioMérieux NucliSENS easyMAG Silica Magnetic Beads

- 175 • bioMérieux NucliSENS easyMAG disposables

- 176 • Biohit pipettor

177

### 178 **Warnings and Precautions**

- 179 • For *In Vitro* Diagnostic Use under Emergency Use Authorization only.

- 180 • Positive results are indicative of the presence of SARS-CoV-2 RNA.

- 181 • Laboratories within the United States and its territories are required to report all positive results to the appropriate public health authorities.

- 183 • The assay has been validated using bioMérieux NucliSENS easyMAG software version 2.0. Please contact Quidel Technical Support prior to modifying or upgrading beyond this version of software.

- 185 • The assay has been validated using Applied Biosystems 7500 Fast Dx software version 1.4. Please contact Quidel Technical Support prior to modifying or upgrading beyond this version of software.

- 187 • The assay has been validated using Applied Biosystems 7500 Standard software version 2.0.6. Please contact Quidel Technical Support prior to modifying or upgrading beyond this version of software.

- 189 • The assay has been validated using Roche LightCycler 480 Instrument II, software version 1.5.0.39. Please contact Quidel Technical Support prior to modifying or upgrading beyond this version of software.

- 191 • The assay has been validated using Qiagen Rotor-Gene Q, software version 2.0.2.4. Please contact Quidel Technical Support prior to modifying or upgrading beyond this version of software.

192

- 193 • The assay has been validated using Bio-Rad CFX96 Touch, software version 3.1. Please contact Quidel  
194 Technical Support prior to modifying or upgrading beyond this version of software.
- 195 • The assay has been validated using Thermo Fisher QuantStudio 7 Pro, software version 2.0. Please contact  
196 Quidel Technical Support prior to modifying or upgrading beyond this version of software.
- 197 • Performance characteristics of this test have been established with the specimen types listed in the  
198 **Intended Use Section** only. The performance of this assay with other specimen types or samples has not  
199 been evaluated.
- 200 • Use of this product should be limited to personnel with sufficient training in PCR and RT-PCR techniques.
- 201 • Treat all specimen/samples as potentially infectious. Follow universal precautions when handling  
202 samples, this kit and its contents.
- 203 • Proper sample collection, storage and transport are essential for correct results.
- 204 • Store assay reagents as indicated on their individual labels.
- 205 • Wear suitable protective clothing, gloves, eye and face protection when using this kit.
- 206 • For accurate results, pipette carefully using only calibrated equipment.
- 207 • Thoroughly clean and disinfect all surfaces with a 10% bleach solution followed by molecular grade water.
- 208 • Use micropipettes with an aerosol barrier or positive displacement tips for all procedures.
- 209 • Avoid microbial and cross contamination of the kit reagents. Follow Good Laboratory Procedures.
- 210 • Do not mix reagents from kits with different lot numbers.
- 211 • Do not use reagents from other manufacturers with this kit.
- 212 • Do not use product after its expiration date.
- 213 • Proper workflow planning is essential to minimize contamination risk. Always plan laboratory workflow in  
214 a uni-directional manner, beginning with pre-amplification and moving through amplification and  
215 detection.
- 216 • Use dedicated supplies and equipment in pre-amplification and amplification areas.
- 217 • Do not allow cross movement of personnel or equipment between areas.
- 218 • Keep amplification supplies separate from pre-amplification supplies at all times.
- 219 • Do not open sample tubes or unseal plates post amplification.
- 220 • Dispose of amplified material carefully and in accordance with local laws and regulations in order to  
221 minimize the risk of amplicon contamination.
- 222 • Do not use supplies dedicated for reagent or sample preparation for processing target nucleic acid.
- 223 • MSDS is available upon request or can be accessed on the product website.

224

### 225 **Storage and Handling of Kit Reagents**

- 226 • Store the unopened kit at 2°C to 8°C until the expiration date listed on the outer kit box.
- 227 • The rehydrated Master Mix may be stored at room temperature (20°C to 25°C) for up to 24 hours. For  
228 longer storage the rehydrated Master Mix should be recapped, sealed with parafilm and stored in an  
229 upright position at ≤-20°C for up to 14 days. Protect the Master Mix from light during storage.

230

### 231 **Indications of Instability or Deterioration of Reagents**

232 Cloudiness of the Rehydration Solution, when within expiration, may indicate deterioration of this reagent.  
233 Contact Quidel Technical Assistance for a replacement.

234

### 235 **Specimen Collection, Storage and Handling**

236 Nasopharyngeal, oropharyngeal, or nasal specimens should be collected, transported, stored, and processed  
237 according to CLSI M41-A.<sup>2</sup> Specimens should be stored at 2°C to 8°C until tested. If specimens cannot be  
238 tested within 72 hours of collection, they should be frozen at -70°C or colder until tested.

239

240 The following viral transport media (M4, M4-RT, M5, M6, MTM and UTM) (1 mL and 3 mL) are compatible  
241 with the Lyra respiratory assays.

242

243 CDC Viral Transport Media ([https://www.cdc.gov/coronavirus/2019-ncov/downloads/Viral-Transport-](https://www.cdc.gov/coronavirus/2019-ncov/downloads/Viral-Transport-Medium.pdf)  
 244 [Medium.pdf](https://www.cdc.gov/coronavirus/2019-ncov/downloads/Viral-Transport-Medium.pdf) ) is compatible with the Lyra SARS-CoV-2 Assay.



245  
 246 **Nucleic Acid Extracts Storage**

247 Eluates from the NucliSENS easyMAG can be stored at room temperature (20°C to 25°C) for 2 hours, at 2°C to  
 248 8°C for 24 hours and 1 month at –20°C to –70°C.

249  
 250 **bioMérieux NucliSENS easyMAG Nucleic Acid Extraction Programming Instructions**



251 Note: A Positive Control (i.e. Lyra SARS-CoV-2, Positive Control #M5153), and a negative process control (i.e.,  
 252 Lyra SARS-CoV-2, Negative Control # M5275) should be included in each extraction run.


253 1. Turn on the instrument and wait for instrument light to appear orange. Then switch on the  
 254 computer/launch NucliSENS easyMAG software. Do not log into software until the light on the instrument  
 255 has turned green.


256 2. Barcode reagents after pressing the ‘Instrument’  and ‘Reagent Inventory’  buttons.


257 3. To enter samples, press the ‘Daily Use’  button, which will default to the ‘Define Request’   
 258 screen. Select the following settings:


- 259 a. Sample ID: Enter the **sample name** using the keyboard.
- 260 b. Matrix: Select **Other** from the drop-down menu
- 261 c. Request: Select **Generic** from the drop-down menu
- 262 d. Volume (mL): Select **0.200** from the drop-down menu
- 263 e. Eluate (µL): Select **50** from the drop-down menu
- 264 f. Type: Primary
- 265 g. Priority: Normal


266 4. Upon pressing the ‘Save’  button, the sample will appear in the ‘Unassigned Sample’ window on  
 267 the left side of the screen. Press the ‘Enter New Extraction Request’   
 268 button, and repeat the process for additional samples. Alternatively multiple samples can be entered by pressing the ‘Auto


269 Create New Extraction Requests’  button.

270 5. Once all samples are created, go to ‘Organize Runs’ by clicking on the  icon near the top of the

271 page. Create a run by pressing the ‘Create Run’  button. Enter a run name or use the default.

272 6. Add samples to the run by using the ‘Auto Fill Run’   
 273 button (auto fills up to 24 samples from the ‘Unassigned Sample list’ on the left hand side of the screen). Alternatively, individual samples can be


274 moved into and out of the run by using the left and right ‘Positioning icons’  after selecting  
 275 the appropriate sample. The sample order within the run can be changed using the ‘Move Extraction

276 Request Up/Down’ buttons .

277 7. Obtain 1 to 3 (for 8 to 24 samples, respectively) sample vessel(s), and add 20 µL of Process Control to each  
 278 sample well used.



279 8. Add 180 µL of each sample to the appropriate well as designated.

280 9. Go to 'Load Run' by pressing the  button near the top of the screen. Insert tips and sample

281 vessel(s) into the instrument

282 10. Enter the barcode(s) of the sample vessel(s)

283 11. Enter the barcode(s) of silica beads to be used

284 12. Close the instrument lid.


285 13. Assign silica beads to samples as follows:

286 a. Click the reagents symbol below number 1 in the picture below. The lot number of the silica

287 beads should appear below the Silica tab at number 2 in the picture below.

288 b. Highlight and select the samples in the run for which beads need to be assigned (in the box

289 containing number 3 in the picture below)

290 c. Click the  positioning icon (below number 4 in the picture below) to assign the silica lot


291 number to the selected samples


292 d. If the bead symbol to the right of number 5 in the picture below is selected, the silica bead lot

293 number should be displayed for each sample



294

295 14. Print work list by touching 'Load Run' icon followed by pressing the 'Print Work List' icon .

296 15. Press the 'Dispense Lysis'  button. The on-board lysis will take approximately 12 minutes to

297 complete.

298 16. For each sample vessel, prepare magnetic particles using the Biohit pipettor and tips for up to eight

299 reactions as follows:

300 a. Using 1 tip and Program 1, aspirate 550 µL nuclease-free water and dispense into a 1.5 mL DNase

301 / RNase free microfuge tube.

302 b. Vortex the magnetic silica. Using 1 tip and Program 1, aspirate 550 µL of magnetic silica, dispense

303 into the water and mix by vortexing.

304 c. Using 1 tip and Program 2, aspirate 1050 µL of the magnetic silica mixture and dispense 25 µL

305 back into the same tube.

306 d. Dispense 125 µL magnetic silica mixture each into 8 wells of an ELISA strip plate. Discard tip.

307 e. After Lysis is complete (NB: the 'Instrument Status' at the bottom of the screen must be 'IDLE!'),

308 using 8 tips and Program 3, aspirate 100 µL of magnetic silica mixture in strip wells, dispense 100

309 µL of magnetic silica mixture in strip wells, and aspirate 100 µL of magnetic silica mixture in strip

310 wells.

- 311 f. Insert tips into liquid within the sample vessels. Aspirate 800  $\mu$ L then dispense 900  $\mu$ L of  
312 magnetic silica mixture back into vessel. Aspirate 1000  $\mu$ L of magnetic silica mixture from vessel  
313 and dispense 1000  $\mu$ L of magnetic silica back into vessel. Repeat aspiration / dispensing of 1000  
314  $\mu$ L two more times.



- 315 17. Close the instrument and press the 'Start' button to begin the run.  
316 18. Upon completion of run, transfer purified nucleic acid to nuclease-free tubes. Eluates from the NucliSENS  
317 easyMAG can be stored at room temperature (20°C to 25°C) for 2 hours, at 2°C to 8°C for 24 hours and 1  
318 month at -20°C to -70°C.  
319

## 320 Assay Procedure

321 Run the following procedures at controlled room temperature of 20°C to 25°C.

322

### 323 Master Mix Rehydration Procedure

- 324 1. Determine the number of specimens extracted to be tested and obtain the correct number of eight-  
325 test lyophilized Master Mix vials for testing.  
326 2. Return unused reagents to the appropriate storage conditions.  
327 3. Open Master Mix carefully to avoid disruption of the pellet.  
328 4. Add 135  $\mu$ L of Rehydration Solution to the Master Mix.  
329 5. Place vial at room temperature for 1 to 2 minutes to allow rehydration of pellet.  
330 6. Gently pipette up and down 2 to 3 times avoiding the formation of bubbles prior to dispensing into  
331 the first plate well or tube.

332 **Note:** The rehydrated Master Mix is sufficient for 8 reactions.

333 **Note:** The rehydrated Master Mix may be stored at room temperature (20°C to 25°C) for up to 24  
334 hours. For longer storage the rehydrated Master Mix should be recapped, sealed with parafilm and  
335 stored in an upright position at  $\leq$ -20°C for up to 14 days. Protect the Master Mix from light during  
336 storage.  
337

### 338 RT-PCR Set-up Procedure

- 339 1. Add 15  $\mu$ L of the rehydrated Master Mix to each plate well or tube.  
340 2. Add 5  $\mu$ L of extracted nucleic acid (specimen with the process control) into the plate well or tube.  
341 Mixing of reagents is not required.  
342 **Note:** Use a new barrier micropipettor tip with each extracted specimen.  
343 3. Seal the plate or tubes.  
344 4. Centrifuge the plate or tubes for a minimum of 15 seconds. Ensure that all liquid is at the bottom of  
345 the plate wells or tubes.  
346 5. Turn on the appropriate thermocycler.  
347 6. Insert plate or tubes into the appropriate thermocycler.  
348

349 **NOTE:** Refer to Appendix for specific programming and testing protocols of each thermocycler.  
350

## 351 Quality Control

352 The Lyra SARS-CoV-2 Assay incorporates several controls to monitor assay performance.  
353

- 354 1. The **Process Control (PRC)** consists of an inactivated and stabilized MS2 Bacteriophage that contains  
355 an RNA genome. It must be used during extraction and amplification in the assay. This control should  
356 be added to each sample aliquot prior to extraction. The PRC serves to monitor inhibitors in the  
357 extracted specimen, assures that adequate amplification has taken place and confirms that the  
358 nucleic acid extraction was sufficient.  
359

- 360 2. The **Positive Control** (containing SARS-CoV-2 Synthetic RNA, Part M5153) must be treated as a  
 361 patient specimen and be included in every extraction and RT-PCR run.  
 362  
 363 3. The **Negative Control** (Part M5275) must be treated as a patient specimen and be included in every  
 364 extraction and PCR run.  
 365  
 366 4. Failure of either the **Positive Control** or the **Negative Control** invalidates the RT-PCR run and results  
 367 should not be reported. The RT-PCR run should be repeated with the extracted controls and  
 368 specimens first. Re-extract and retest another aliquot of the controls and the specimens or obtain  
 369 new samples and retest if the controls fail again.  
 370

**Table 3. Expected Results from Controls (Applied Biosystems 7500 Fast Dx, Applied Biosystems 7500 Standard, Bio-Rad CFX96 Touch, Qiagen Rotor-Gene Q, or Thermo Fisher QuantStudio 7 Pro)**

Control Type/ Name	Used to Monitor	SARS-CoV-2	Expected Ct Values	PRC	Expected Ct Values
<b>Positive Control</b>	Substantial reagent failure including primer and probe integrity	+	5.0 ≤ Ct ≤ 30.0	+/-	NA <sup>1</sup>
<b>Negative Control</b>	Reagent and/or environmental contamination	-	None detected	+	5.0 ≤ Ct ≤ 30.0

371

**Expected Results from Controls (Roche LightCycler 480)**

Control Type/ Name	Used to Monitor	SARS-CoV-2	Expected Ct Values	PRC	Expected Ct Values
<b>Positive Control</b>	Substantial reagent failure including primer and probe integrity	+	5.0 ≤ Ct ≤ 40.0	+/-	NA <sup>1</sup>
<b>Negative Control</b>	Reagent and/or environmental contamination	-	None detected	+	5.0 ≤ Ct ≤ 40.0

372 <sup>1</sup>No Ct value is required for the Process Control to make a positive call.

373

374 **Interpretation of Results from Patient Specimens**

375

**Table 4. Interpretation of the Lyra SARS-CoV-2 Assay Results on the Applied Biosystems 7500 Fast Dx, Applied Biosystems 7500 Standard, Bio-Rad CFX96 Touch, Qiagen Rotor-Gene Q, or Thermo Fisher QuantStudio 7 Pro**

Assay Result	Detector: SARS-CoV-2	Detector: Process Control	Interpretation of Results	Notes and Special Guidance
Negative	No Ct detected	5.0 ≤ Ct ≤ 30.0	No SARS-CoV-2 viral RNA detected; PRC Detected.	
SARS-CoV-2 Positive	5.0 ≤ Ct ≤ 30.0	NA <sup>1</sup>	SARS-CoV-2 Virus viral RNA detected.	

Invalid	No Ct detected	No Ct detected	No SARS-CoV-2 viral RNA and no PRC RNA detected.	Invalid test. Retest the same purified sample. If the test is also invalid, re-extract and retest another aliquot of the same specimen or obtain a new specimen and retest.
<b>Interpretation of the Lyra SARS-CoV-2 Assay Results on the Roche LightCycler 480</b>				
<b>Assay Result</b>	<b>Detector: SARS-CoV-2</b>	<b>Detector: Process Control</b>	<b>Interpretation of Results</b>	<b>Notes and Special Guidance</b>
Negative	No Ct detected	5.0 ≤ Ct ≤ 40.0	No SARS-CoV-2 viral RNA detected; PRC Detected.	
SARS-CoV-2 Positive	5.0 ≤ Ct ≤ 40.0	NA <sup>1</sup>	SARS-CoV-2 Virus viral RNA detected.	
Invalid	No Ct detected	No Ct detected	No SARS-CoV-2 viral RNA and no PRC RNA detected.	Invalid test. Retest the same purified sample. If the test is also invalid, re-extract and retest another aliquot of the same specimen or obtain a new specimen and retest.

<sup>1</sup> No Ct value is required for the Process Control to make a positive call.

376  
377  
378  
379  
380  
381  
382  
383  
384  
385  
386  
387  
388  
389

**CLINICAL PERFORMANCE**

The clinical performance of the Lyra SARS-CoV-2 Assay was evaluated using two different studies:

- A study using two hundred sixty-five fresh or frozen nasopharyngeal swab specimens collected in UTM (36 and 229, respectively) from patients located in the USA.
- A fully contrived positive specimen study using nasopharyngeal swab specimens.

All two hundred sixty-five specimens were negative for SARS-CoV-2 when extracted with the NucliSENS easyMAG system and tested by the Lyra SARS-CoV-2 Assay.

One hundred twenty-four specimens included in this study were positive for other circulating respiratory viruses as identified by FDA-cleared assays:

<b>Circulating Virus</b>	<b># of positive specimens</b>
Influenza A	30
RSV	4
Coronavirus Seasonal (not identified)	10
Coronavirus 229e	20
Coronavirus OC43	20
Coronavirus NL63	20
Coronavirus HKU1	20

390  
391  
392  
393  
394  
395  
396

Viral RNA was obtained from BEI Resources Repository for use in the contrived clinical study. The genomic RNA was extracted from a preparation of cell lysate and supernatant from Cercopithecus aethiops kidney epithelial (Vero E6, ATCC® CRL-1586™) cells infected with SARS-related coronavirus 2 (SARS-CoV-2), isolate USA-WA1/2020, using Qiagen QIAamp® Viral RNA Mini Kit (Qiagen 52904). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA. The Genome Copy Number was established using Bio-Rad QX200 Droplet Digital PCR (ddPCR™) System.

397

398

Ninety-two positive contrived samples were created by spiking ninety-two individual clinical samples determined to be negative for SARS-CoV-2 by the Lyra SARS-CoV-2 Assay prior to spiking with one of three concentrations of genomic SARS-CoV-2 RNA. Forty-four specimens were spiked with 1x LoD (8.00E-01 cp/μL) of RNA. Twenty-four additional specimens were spiked with 3x LoD (2.40E00 cp/μL) of RNA. Twenty-four additional specimens were spiked with 5x LoD (4.00E00 cp/μL) of RNA. All samples were extracted and tested according to the Lyra SARS-CoV-2 Assay Package Insert.

404

405

All ninety-two contrived samples were positive in the Lyra SARS-CoV-2 Assay. The results for the contrived positive specimens are shown in the table below:

406

407

Sample RNA Concentration	# Positives/# Tested	Mean SARS-CoV-2 Ct	%CV
unspiked	0/92	NA	NA
1 .0x LoD	44/44	26.9	5.7
3x LoD	24/24	22.8	3.4
5x LoD	24/24	22.4	3.0

408

Performance against the expected results are:

409

Positive Percent Agreement 92/92 = 100% (95% CI: 95.99%-100%)

410

Negative Percent Agreement 92/92 = 100% (95% CI: 95.99%-100%)

411

412

**ANALYTICAL PERFORMANCE**

413

**Level of Detection**

414

The Limit of Detection of the Lyra SARS-CoV-2 Assay utilized limiting dilutions of genomic SARS-CoV-2 RNA in negative nasopharyngeal matrix. Each dilution was extracted using the NucliSENS easyMAG System and tested on Applied Biosystems 7500 Fast Dx, Applied Biosystems 7500 Standard, Roche LightCycler 480, Qiagen Rotor-Gene Q, Bio-Rad CFX96 Touch, or Thermo Fisher QuantStudio 7 Pro. Analytical sensitivity (LoD) is defined as the lowest concentration at which at least 95% of all replicates tested positive.

418

419

420

The genomic RNA was extracted from a preparation of cell lysate and supernatant from Cercopithecus aethiops kidney epithelial (Vero E6, ATCC® CRL-1586™) cells infected with SARS-related coronavirus 2 (SARS-CoV-2), isolate USA-WA1/2020, using Qiagen QIAamp Viral RNA Mini Kit (Qiagen 52904). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA. The Genome Copy Number was established using Bio-Rad QX200 Droplet Digital PCR (ddPCR) System.

424

425

426

This study established the LoD for the Lyra SARS-CoV-2 Assay as 8.00E-01 genomic RNA copies/μL, subsequently confirmed by testing 20 replicates.

427

428

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	1	23.95	18.54	Positive
	2	26.59	18.28	Positive
	3	26.19	18.32	Positive
	4	25.13	18.41	Positive
	5	24.88	18.74	Positive
	6	24.84	19.18	Positive
	7	25.51	18.82	Positive
	8	25.20	18.58	Positive

**Table 6. LoD in Nasopharyngeal specimens with Applied Biosystems 7500 Fast Dx**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
	9	24.69	18.71	Positive
	10	24.57	18.67	Positive
	11	23.86	18.75	Positive
	12	24.58	18.91	Positive
	13	25.19	19.03	Positive
	14	25.84	19.05	Positive
	15	26.58	19.10	Positive
	16	26.72	19.15	Positive
	17	24.16	19.06	Positive
	18	25.15	18.91	Positive
	19	25.51	19.05	Positive
	20	24.41	19.07	Positive

429

**Table 7. LoD in Oropharyngeal specimens with Applied Biosystems 7500 Fast Dx**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	1	27.26	19.38	Positive
	2	28.99	19.22	Positive
	3	27.3	19.51	Positive
	4	26.09	19.27	Positive
	5	26.88	19.61	Positive
	6	26.02	19.19	Positive
	7	26.37	19.21	Positive
	8	25.01	19.30	Positive
	9	25.14	19.06	Positive
	10	26.21	19.03	Positive
	11	27.79	19.27	Positive
	12	28.83	19.12	Positive
	13	28.83	19.19	Positive
	14	26.81	19.50	Positive
	15	25.1	19.30	Positive
	16	26.2	19.11	Positive
	17	26.74	19.00	Positive
	18	25.28	19.13	Positive
	19	26.27	19.31	Positive
	20	26.37	19.24	Positive

430

**Table 8. LoD in Nasopharyngeal specimens with Applied Biosystems 7500 Standard**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	1	26.63	19.26	Positive
	2	29.15	19.28	Positive
	3	25.67	19.69	Positive
	4	25.53	20.07	Positive
	5	26.15	20.50	Positive
	6	26.71	20.50	Positive

**Table 8. LoD in Nasopharyngeal specimens with Applied Biosystems 7500 Standard**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
	7	26.11	19.14	Positive
	8	26.94	19.18	Positive
	9	25.62	18.64	Positive
	10	25.80	18.80	Positive
	11	26.76	19.15	Positive
	12	26.15	19.63	Positive
	13	27.42	19.44	Positive
	14	27.51	19.99	Positive
	15	26.07	19.9	Positive
	16	25.92	18.81	Positive
	17	27.95	20.02	Positive
	18	27.71	19.27	Positive
	19	26.51	18.86	Positive
	20	Undetermined	19.11	Negative

431

**Table 9. LoD in Nasopharyngeal specimens with Roche LightCycler 480\***

Concentration	Replicate	SARS-CoV-2 Ct*	PRC Ct*	Interpretation
8.00E-01 genomic RNA copies/μL	1	32.91	31.73	Positive
	2	34.54	32.9	Positive
	3	34.83	32.25	Positive
	4	34.94	31.7	Positive
	5	33.81	32.14	Positive
	6	34.36	32.37	Positive
	7	33.90	32.10	Positive
	8	33.83	32.80	Positive
	9	33.8	31.86	Positive
	10	34.28	32.27	Positive
	11	33.63	32.81	Positive
	12	33.72	32.45	Positive
	13	34.86	33.17	Positive
	14	34.57	32.64	Positive
	15	34.48	32.92	Positive
	16	33.61	32.82	Positive
	17	33.87	33.34	Positive
	18	34.44	33.36	Positive
	19	34.22	32.55	Positive
	20	33.77	32.97	Positive

432

\* Results include 10 cycles not captured by the other instruments

433

**Table 10. LoD in Nasopharyngeal specimens with Qiagen Rotor-Gene Q**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	1	24.01	19.08	Positive
	2	24.04	19.36	Positive
	3	24.85	19.44	Positive

**Table 10. LoD in Nasopharyngeal specimens with Qiagen Rotor-Gene Q**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
	4	23.23	19.13	Positive
	5	24.39	19.07	Positive
	6	23.89	18.94	Positive
	7	23.78	18.80	Positive
	8	24.82	18.86	Positive
	9	23.87	18.83	Positive
	10	24.05	18.90	Positive
	11	23.28	18.84	Positive
	12	24.36	18.71	Positive
	13	23.85	18.87	Positive
	14	23.54	18.88	Positive
	15	24.84	19.20	Positive
	16	23.63	19.01	Positive
	17	24.18	18.97	Positive
	18	23.47	19.01	Positive
	19	23.58	18.94	Positive
	20	23.89	19.02	Positive

434

**Table 11. LoD in Nasopharyngeal specimens with Bio-Rad CFX96 Touch**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	1	27.19	21.25	Positive
	2	25.57	21.35	Positive
	3	25.80	22.68	Positive
	4	27.93	21.3	Positive
	5	29.03	21.09	Positive
	6	25.79	21.45	Positive
	7	25.65	21.19	Positive
	8	26.26	21.16	Positive
	9	29.46	21.41	Positive
	10	25.09	21.45	Positive
	11	25.68	21.36	Positive
	12	28.51	21.49	Positive
	13	25.5	21.97	Positive
	14	26.81	21.36	Positive
	15	26.17	21.1	Positive
	16	25.04	21.91	Positive
	17	25.47	22.08	Positive
	18	25.54	21.26	Positive
	19	25.77	22.29	Positive
	20	25.59	22.16	Positive

435

**Table 12. LoD in Nasopharyngeal specimens with Thermo Fisher QuantStudio 7 Pro**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
	1	24.25	20.21	Positive



**Table 12. LoD in Nasopharyngeal specimens with Thermo Fisher QuantStudio 7 Pro**

Concentration	Replicate	SARS-CoV-2 Ct	PRC Ct	Interpretation
8.00E-01 genomic RNA copies/μL	2	26.7	20.9	Positive
	3	27.14	20.6	Positive
	4	27.28	20.81	Positive
	5	29.60	20.78	Positive
	6	26.99	20.65	Positive
	7	28.75	20.82	Positive
	8	27.63	20.76	Positive
	9	29.80	20.65	Positive
	10	26.60	20.55	Positive
	11	27.23	20.54	Positive
	12	29.81	20.73	Positive
	13	26.59	20.88	Positive
	14	27.23	20.87	Positive
	15	26.63	20.62	Positive
	16	26.07	20.84	Positive
	17	25.14	20.81	Positive
	18	27.34	20.6	Positive
	19	29.22	20.67	Positive
	20	26.37	20.38	Positive

436

437

**Analytical Reactivity (Inclusivity)**

438

The inclusivity of the Lyra SARS-CoV-2 Assay was established by testing Genomic RNA from the SARS-related coronavirus 2 (SARS-CoV-2), isolate USA-WA1/2020, via *in-silico* analysis. The *in-silico* analysis demonstrated the Lyra SARS-CoV-2 primers are 100% conserved to 257 SARS-CoV-2 sequences available from NCBI and GISAID as of March 5, 2020.

440

441

442

**Analytical Specificity (Cross-Reactivity)**

443

The Analytical Specificity of the assay was established by both direct testing of organisms in the assay (“wet” testing) and *in silico* analysis. The wet testing used 25 micro-organisms, in high concentrations, identified by the FDA as high priority for evaluation due to the reasonable likelihood they may be present in upper respiratory samples. All micro-organisms were undetectable with the Lyra SARS-CoV-2 Assay when wet tested as shown below.

444

445

446

447

448

449

**Table 13. Cross-reactivity test results**

Virus/Bacteria/Parasite	Strain	Source/ Sample type	Concentration	Results
Adenovirus	Type 1	Isolate	1 x 10 <sup>7.53</sup> U/mL	Neg, Neg, Neg
Coronavirus	229e	Isolate	1 x 10 <sup>6.10</sup> U/mL	Neg, Neg, Neg
Coronavirus	OC43	Isolate	9.55 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Neg, Neg, Neg
Coronavirus	NL63	Isolate	1 x 10 <sup>4.67</sup> U/mL	Neg, Neg, Neg
MERS-CoV (heat-inactivated)	Florida/USA- 2_Saudia Arabia_2014	Isolate	4.17 x 10 <sup>5</sup> TCID <sub>50</sub> /mL	Neg, Neg, Neg
SARS -1	2003-00592	Inactivated virus	Not available	Neg, Neg, Neg
<i>Mycoplasma pneumoniae</i>	M129	Isolate	3 x 10 <sup>7</sup> CCU/mL	Neg, Neg, Neg

Virus/Bacteria/Parasite	Strain	Source/ Sample type	Concentration	Results
<i>Streptococcus pyogenes</i>	Z018	Isolate	3.8 x 10 <sup>9</sup> cfu/mL	Neg, Neg, Neg
Influenza A H3N2	Brisbane/10/07	Isolate	1 x 10 <sup>5.07</sup> U/mL	Neg, Neg, Neg
Influenza A H1N1	New Caledonia/20/99	Isolate	1 x 10 <sup>6.66</sup> U/mL	Neg, Neg, Neg
Influenza B	Brisbane/33/08	Isolate	1 x 10 <sup>5.15</sup> U/mL	Neg, Neg, Neg
Parainfluenza	Type 1	Isolate	1 x 10 <sup>8.01</sup> U/mL	Neg, Neg, Neg
Parainfluenza	Type 2	Isolate	1 x 10 <sup>6.34</sup> U/mL	Neg, Neg, Neg
Parainfluenza	Type 3	Isolate	8.51 x 10 <sup>7</sup> TCID50/mL	Neg, Neg, Neg
Parainfluenza	Type 4b	Isolate	1 x 10 <sup>7.53</sup> U/mL	Neg, Neg, Neg
Enterovirus	Type 68	Isolate	1 x 10 <sup>6.5</sup> U/mL	Neg, Neg, Neg
Human Metapneumovirus	A1 (IA10-s003)	Isolate	1 x 10 <sup>5.55</sup> U/mL	Neg, Neg, Neg
Respiratory Syncytial Virus	Type A (3/2015 Isolate #3)	Isolate	1 x 10 <sup>5.62</sup> U/mL	Neg, Neg, Neg
Human Rhinovirus	N/A	Inactivated virus	Not available	Neg, Neg, Neg
<i>Chlamydomphila pneumoniae</i>	AR-39	Isolate	2.9 x 10 <sup>7</sup> IFU/mL	Neg, Neg, Neg
<i>Haemophilus influenzae</i>	Type b; Eagan	Isolate	7.87 x 10 <sup>8</sup> cfu/mL	Neg, Neg, Neg
<i>Legionella pneumophila</i>	Philadelphia	Isolate	6.82 x 10 <sup>9</sup> cfu/mL	Neg, Neg, Neg
<i>Streptococcus pneumoniae</i>	Z022; 19f	Isolate	2.26 x 10 <sup>9</sup> cfu/mL	Neg, Neg, Neg
<i>Bordetella pertussis</i>		Isolate		Neg, Neg, Neg
<i>Pneumocystis jirovecii</i> -S. <i>cerevisiae</i> Recombinant	W303-Pji	Isolate	1.56 x 10 <sup>8</sup> cfu/mL	Neg, Neg, Neg
Negative Nasopharyngeal Matrix	MTM	N/A	N/A	Neg, Neg, Neg
Negative Nasopharyngeal Matrix	MTM	N/A	N/A	Neg, Neg, Neg
Negative Nasal Matrix	CDC Viral Transport	N/A	N/A	Neg, Neg, Neg
Negative Oropharyngeal Matrix	CDC Viral Transport	N/A	N/A	Neg, Neg, Neg

450

451 The *in silico* analysis focused on 32 micro-organisms identified by the FDA as high priority for  
452 assessment due to their potential presence in upper respiratory samples.

453

Organism	Total # Sequences	# Complete Genomes	# WGS Strains
Adenovirus	532	532	0
Coronavirus (Seasonal)	288	288	0
Enterovirus <sup>B</sup>	2708	2674	34
Influenza A Virus <sup>A B</sup>	172455	21444 (+39 A/Mexico/4108/2009)	108
Influenza B Virus <sup>A B</sup>	53952	6755 (+16 B/Florida/4/2006)	0
Influenza C Virus <sup>B</sup>	2205	N/A	N/A
Human Metapneumovirus	145	145	0
Human Parainfluenza Virus 1-4	439	439	0
Human Parechovirus	124	124	0

<b>Organism</b>	<b>Total # Sequences</b>	<b># Complete Genomes</b>	<b># WGS Strains</b>
Human Respiratory Syncytial Virus <sup>B</sup>	1275	1275	0
Rhinovirus	214	214	0
SARS-1	236 <sup>C</sup>	232 (+4 pp1ab sequences)	0
<i>Bacillus anthracis</i>	4152	69	86
<i>Candida albicans</i>	1541	59	34
<i>Chlamydia pneumoniae</i>	466	5	20
<i>Chlamydia psittaci</i>	11179	23	45
<i>Corynebacterium diphtheriae</i>	20797	17	194
<i>Coxiella burnetii</i>	419	28	3
<i>Haemophilus influenzae</i>	45267	61	692
Legionella <sup>B</sup>	4843	98	65
Leptospira <sup>B</sup>	64456	133	266
<i>Moraxella catarrhalis</i> <sup>B</sup>	8333	11	184
<i>Mycobacterium tuberculosis</i>	194	194	0
<i>Mycoplasma pneumoniae</i>	808	51	45
<i>Neisseria elongata</i> & <i>N. meningitidis</i> <sup>B</sup>	312050	116	1318
<i>Pneumocystis jirovecii</i>	487	15	3
<i>Pseudomonas aeruginosa</i>	195	195	0
<i>Staphylococcus aureus</i>	634	634	0
<i>Staphylococcus epidermidis</i> <sup>B</sup>	61880	23	508
<i>Streptococcus pneumoniae</i> <sup>B</sup>	1633369	107	8526
<i>Streptococcus pyogenes</i> <sup>B</sup>	46153	201	1733
<i>Streptococcus salivarius</i> <sup>B</sup>	9417	18	48
<sup>A</sup> Genome counts for Influenza A and Influenza B were attained for strains that included all 8 segments, except for A/Mexico/4108/2009(H1N1) and B/Florida/4/2006; all available gene sequences were included.			
<sup>B</sup> For BLAST, 'Max Target Seqs' was set to 5000. See Table 2.			
<sup>C</sup> 4 polyprotein cds sequences were also included.			

454

455

456

457

458

459

460

461

The in-silico analysis demonstrated < 80% homology with all organisms except for the following: three Enterovirus sequences are 80.9% conserved to the reverse primer, however, the forward primer is only 76% conserved and the probe alignment had an overall homology of 56%. The SARS-1 sequences are ≥80% conserved to both primers, however, the last base on the 3' ends of both primers are not conserved. The wet testing of the only available SARS-1 strain was non-detectable.

**Limitations**

462

463

464

465

466

467

- Negative results do not preclude infection with SARS-CoV-2 and should not be the sole basis of a patient treatment decision.
- This test is intended to be used for the detection of SARS-CoV-2 RNA in nasopharyngeal and oropharyngeal swab samples. Testing of other sample types may result in inaccurate results.
- Nasal swabs and mid-turbinate nasal swabs are considered acceptable specimen types for use with the Lyra SARS-CoV-2 Assay but performance with these specimen types has not been established. Testing of

468 nasal and mid-turbinate nasal swabs (self-collected under supervision of or collected by a healthcare  
469 provider) is limited to patients with symptoms of COVID-19. Please refer to FDA's FAQs on Diagnostic  
470 Testing for SARS-CoV-2 for additional information.

- 471 • Improper collection, storage or transport of specimens may lead to false negative results.
- 472 • Inhibitors present in the sample and/or errors in following the assay procedure may lead to false negative  
473 results.
- 474 • A trained health care professional should interpret assay results in conjunction with the patient's medical  
475 history, clinical signs and symptoms, and the results of other diagnostic tests.
- 476 • Analyte targets (viral sequences) may persist *in vivo*, independent of virus viability. Detection of analyte  
477 target(s) does not imply that the corresponding virus(es) are infectious, nor are the causative agents for  
478 clinical symptoms.
- 479 • There is a risk of false positive values resulting from cross-contamination by target organisms, their  
480 nucleic acids or amplified product, or from non-specific signals in the assay.
- 481 • There is a risk of false negative values due to the presence of sequence variants in the viral targets of the  
482 assay.
- 483 • The assay performance was not established in immunocompromised patients.

484

#### 485 **Conditions of Authorization for the Labs**

486 The Lyra SARS-CoV-2 Assay Letter of Authorization, along with the authorized Fact Sheet for Healthcare  
487 Providers, the authorized Fact Sheet for Patients, and authorized labeling are available on the FDA website:  
488 <https://www.fda.gov/MedicalDevices/Safety/EmergencySituations/ucm161496.htm>.

489 However, to assist clinical laboratories using the Lyra SARS-CoV-2 Assay, the relevant Conditions of  
490 Authorization are listed below.

- 491 • Authorized laboratories<sup>1</sup> using the Lyra SARS-CoV-2 Assay will include with result reports of the Lyra  
492 SARS-CoV-2 Assay test, all authorized Fact Sheets. Under exigent circumstances, other appropriate  
493 methods for disseminating these Fact Sheets may be used, which may include mass media.
- 494 • Authorized laboratories using the Lyra SARS-CoV-2 Assay will perform the Lyra SARS-CoV-2 Assay as  
495 outlined in the Lyra SARS-CoV-2 Assay Instructions for Use. Deviations from the authorized  
496 procedures, including the authorized instruments, authorized extraction methods, authorized clinical  
497 specimen types, authorized control materials, authorized other ancillary reagents and authorized  
498 materials required to perform the Lyra SARS-CoV-2 Assay are not permitted.
- 499 • Authorized laboratories that receive the Lyra SARS-CoV-2 Assay must notify the relevant public health  
500 authorities of their intent to run the test prior to initiating testing.
- 501 • Authorized laboratories using the Lyra SARS-CoV-2 Assay will have a process in place for reporting  
502 test results to healthcare providers and relevant public health authorities, as appropriate.
- 503 • Authorized laboratories will collect information on the performance of the test and report to  
504 DMD/OHT7-OIR/OPEQ/CDRH (via email: [CDRH-EUA-Reporting@fda.hhs.gov](mailto:CDRH-EUA-Reporting@fda.hhs.gov)) and Quidel  
505 (QDL.COV2.test.event.report@quidel.com) any suspected occurrence of false positive or false  
506 negative results and significant deviations from the established performance characteristics of the  
507 test of which they become aware.
- 508 • All laboratory personnel using the test must be appropriately trained in RT-PCR techniques and use  
509 appropriate laboratory and personal protective equipment when handling this kit, and use the test in  
510 accordance with the authorized labeling.
- 511 • Quidel, its authorized distributor(s) and authorized laboratories using the Lyra SARS-CoV-2 Assay will  
512 ensure that any records associated with this EUA are maintained until otherwise notified by FDA.  
513 Such records will be made available to FDA for inspection upon request.
- 514 •

514

515 <sup>1</sup>For ease of reference, the letter of authorization refers to, "United States (U. S.) laboratories certified under the  
516 Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, to perform high complexity tests" as  
517 "authorized laboratories."  
518

**519 Customer and Technical Assistance**

520 To place an order or for technical support, please contact a Quidel Representative at 800.874.1517 (toll-free in  
521 the U.S.) or 858.552.1100 (outside of U.S.), Monday through Friday, between 8:00 a.m. and 5:00 p.m., Eastern  
522 Time. Orders may also be placed by fax at 740.592.9820. For e-mail support contact:

523 [customerservice@quidel.com](mailto:customerservice@quidel.com) or [technicalsupport@quidel.com](mailto:technicalsupport@quidel.com). For services outside the U.S., please contact  
524 your local distributor. Additional information about Quidel, our products, and our distributors can be found on  
525 our website [quidel.com](http://quidel.com).

526

**527 Intellectual Property**

528 Dye compounds in this product are sold under license from BioSearch Technologies, Inc., and protected by  
529 U.S. and world-wide patents either issued or under application.

530

531 Quasar and FAM are trademarks of Biosearch. NucliSENS, easyMAG and EMAG are trademarks of bioMérieux.  
532 TaqMan and LightCycler are trademarks of Roche. Applied Biosystems and QuantStudio are trademarks of  
533 Thermo Fisher Scientific. Rotor-Gene Q and QIAamp are trademarks of Qiagen. CFX96 Touch and ddPCR are  
534 trademarks of Bio-Rad Laboratories.

535

**536 References**

- 537 1. Baker, S., Frias, L., and Bendix, A. Coronavirus live updates: More than 92,000 people have been infected  
538 and at least 3,100 have died. The US has reported 6 deaths. Here's everything we know. Business Insider.  
539 March 03, 2020.
- 540 2. Clinical and Laboratory Standards Institute. Viral Culture; Approved Guidelines. CLSI document M41-A  
541 [ISBN 1562386239] Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne,  
542 Pennsylvania 19087-1898, USA 2006.
- 543 3. Lauer, S.A., et. al. The incubation period of Coronavirus disease 2019 (COVID-19) from publicly reported  
544 confirmed cases: estimation and application, Ann Intern Med. 2020
- 545 4. [www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html](http://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html)

546

547

548

**APPENDIX**

549

**Applied Biosystems 7500 Fast Dx Programming Instructions**

550

Refer to User Manual Part Number 4406991 for additional information.

551

1. Launch the 7500 Fast Dx software package.

552

2. The **Quick Startup document** dialog window will open. Select the **Create New Document** button to start the **New Document Wizard**. Follow each step to initiate the Lyra™ SARS-CoV-2 Assay protocol.

553

554

a. Define Document: Most of the following should be the default setting. If not, change accordingly.

555

i. Confirm or enter the following information.

<b>Assay:</b>	Standard Curve (Absolute Quantitation)
<b>Container:</b>	96-Well Clear
<b>Template:</b>	Blank Document
<b>Run Mode:</b>	Fast 7500
<b>Operator:</b>	<i>your operator name</i>
<b>Comments:</b>	SDS v1.4
<b>Plate Name:</b>	'Lyra SARS-CoV-2 Assay'

556

ii. Select the **Next** button.

557

558

b. Select Detectors: New detectors for SARS-CoV-2 and the process control (PRC) must be added. For each target, select the **New Detector** button to open the **New Detector** pop-up window. Alternatively, use the **Create Another** button from within the **New Detector** pop-up window for the last two detectors.

559

560

561

562

563

i. Enter the following information for each detector.

Name	Reporter Dye	Quencher Dye	Color
SARS-CoV-2	FAM	(none)	(Select)
PRC	Quasar 670	(none)	(Select)

564

565

ii. Select a unique color to represent each detector.

566

iii. Highlight the new detectors and add to the **Detectors in Document** column using the **Add** button.

567

568

iv. Select **(none)** from the **Passive Reference** drop-down menu.

569

v. Select the **Next** button.

570

vi. Select the **Finish** button without setting any wells.

571

c. The wizard will close and the software will open, starting with the **Setup** tab. This will show the sample plate that was set up during the quick start. For the initial set up, nothing needs to be changed here.

572

573

574

d. Defining the Thermocycler Protocol: Select the **Instrument** tab to set up the Lyra™ SARS-CoV-2 Assay RT-PCR cycling times and temperatures. Under **Thermal Profile** there should be a default 2-stage protocol. Each stage will have 3 user-editable text boxes. The top box value represents the number of reps or cycles for that stage. The middle box value represents the temperature (°C) and the lowest box value represents the time (minutes: seconds).

575

576

577

578

579

i. Make the following changes to the default **Thermal Cycler Protocol**:

580

1. Stage 1

581

a. Reps: 1

582

b. Temp: 55

583

c. Time: 5:00

584

2. Select the bar between Stage 1 and Stage 2. Select the **Add Hold** button to add another stage.

585

586

3. Stage 2

587

a. Reps: 1

588

b. Temp: 60

589

- 590 c. Time: 5:00
- 591 4. Select the bar between Stage 2 and Stage 3. Select the **Add Hold** button to add
- 592 another stage.
- 593 5. Stage 3
- 594 a. Reps: 1
- 595 b. Temp: 65
- 596 c. Time: 5:00
- 597 6. Stage 4 (2-Step Dissociation Stage)
- 598 a. Reps: 10
- 599 b. Step 1
- 600 i. Temp: 92
- 601 ii. Time: 0:05
- 602 c. Step 2
- 603 i. Temp: 57
- 604 ii. Time: 0:40
- 605 7. Select the bar to the right of Stage 4. Select the **Add Cycle** button to add another
- 606 stage.
- 607 8. Stage 5 (2-Step Dissociation Stage)
- 608 a. Reps: 30
- 609 b. Step 1
- 610 i. Temp: 92
- 611 ii. Time: 0:05
- 612 c. Step 2
- 613 i. Temp: 57
- 614 ii. Time: 0:40
- 615 9. If a wrong stage is added the stage can be removed by pressing the **Delete** button
- 616 after highlighting the stage between the vertical lines
- 617 ii. Under **Settings** enter the following:

Sample Volume (µL):	20 (default)
Run Mode:	7500 Fast (default)
Data Collection:	Stage 5, Step 2(57.0 @ 0:40)
<b>NOTE: Do not check the check box next to 'Expert Mode'.</b>	

- 618
- 619 e. Set threshold for each analyte.
- 620 i. Select the **Results** tab.
- 621 ii. Select the **Amplification Plot** tab.
- 622 iii. Select SARS-CoV-2 from the Detector tab in the top right corner.
- 623 iv. In the **Analysis Settings** block, set the **Threshold** to **7.5e+004**.
- 624 v. Select the **Auto Baseline** radio button.
- 625 vi. Select PRC from the Detector tab in the top right corner.
- 626 vii. In the **Analysis Settings** block, set the **Threshold** to **1.0e+004**.
- 627 viii. Select the **Auto Baseline** radio button.
- 628
- 629 f. Save the new protocol as a template for future use.
- 630 i. At the top of the screen select **File** and then **Save As**.
- 631 ii. **Save In:** D:\Applied Biosystems\7500 Fast System\Templates\  
632 **File name:** 'Lyra SARS-CoV-2'
- 633 iii. **Save as type:** 'SDS Templates (\*.sdt)'
- 634 iv. **Save as type:** 'SDS Templates (\*.sdt)'
- 635 g. Exit the software.

### Applied Biosystems 7500 Fast Dx Thermocycler Test Procedure

- 636 1. Launch the 7500 Fast Dx software v1.4 package.
- 637

638  
639  
640

2. The **Quick Startup document** dialog window will open.
3. Click on **Create a new document**.
4. Most of the following should be the default setting. If not, change accordingly.

<b>Assay:</b>	Standard Curve (Absolute Quantitation)
<b>Container:</b>	96-Well Clear
<b>Template:</b>	Lyra SARS-CoV-2
<b>Run Mode:</b>	Fast 7500
<b>Operator:</b>	<i>your operator name</i>
<b>Comments:</b>	SDS v1.4
<b>Plate Name:</b>	YYMMDD- Lyra SARS-CoV-2

641  
642  
643  
644  
645  
646  
647  
648  
649  
650  
651  
652  
653  
654  
655  
656  
657  
658  
659  
660  
661  
662  
663

5. Set Up Sample Plate
  - a. Under the **Setup** and **Plate** tabs the plate setup will appear.
  - b. Select all wells that will contain sample, right-click and select the **Well Inspector** from the drop-down menu. When the **Well Inspector** pop-up window opens, select the detectors for SARS-CoV-2 and PRC.
  - c. Use the **Well Inspector** to enter the sample names. Patient IDs can be entered in the Well Inspector window. However, it is recommended that this is done prior to re-suspending the lyophilized master mix, post run or using the import function to minimize the time the PCR reactions will sit at room temperature prior to starting the run.
  - d. Save the run as **YYMMDD- Lyra SARS-CoV-2.sds**.
  - e. A window will open asking for the "Reason for change of entry". Enter "**Setup**" and any other comments relevant to the run.
6. Starting the PCR
  - a. Select the **Instrument** tab.
  - b. Insert the 96 well PCR plate into the machine.
  - c. Under **Instrument Control**, select the **Start** button to initiate the run.
7. Post PCR
 

**IMPORTANT:** When the run is finished press OK.

  - a. Analyze the data by pressing the "**Analyze**" button in the top menu and save the file.
  - b. Save the file by pressing **Save Document** in the task bar. A window will open asking for the "Reason for change of entry".
  - c. Enter "**Data analysis post run**" and any other comments relevant to the run.

**Applied Biosystems 7500 Standard Programming Instructions**

Refer to User Manual Part Number 4387783 rev C for additional information.

666  
667  
668  
669  
670  
671  
672  
673  
674  
675  
676  
677

1. Launch the 7500 software package.
2. Select the **Advanced Setup** button to open Setup and Experiment Properties. Follow each step to initiate the Lyra SARS-CoV-2 protocol.
  - a. Experiment Name: Enter Experiment Name as SARS-CoV-2. Leave the Barcode, User Name, and Comments fields blank
  - b. Define Experiment Setup: Select 7500 (96 Wells), Quantitation- Standard Curve, TaqMan® Reagents, and Standard (~2 hours to complete a run)
3. In the upper left menu select **Plate Setup**
  - a. Define Targets: New detectors for SARS-CoV-2, and the process control (PRC) must be added.
    - i. Enter the following information for each detector.

Name	Reporter Dye	Quencher Dye	Color
SARS-CoV-2	FAM	(none)	(Select)
PRC	Quasar 670	(none)	(Select)

678



- 
- 679                   ii. Select **Add New Target** button for each target.
- 680                   iii. From each drop down menu select reporter, quencher, and color
- 681                   iv. Select a unique color to represent each detector
- 682           b. Assign Targets and Samples: Under this tab in the bottom left corner, select **none** as the
- 683           Passive Reference.
- 684   4. Select **Run Method** from the upper left menu
- 685           a. Set the **Reaction Volume** per Well to 20 µL under the **Graphical** or **Tabular View**
- 686           b. Define the Thermocycler Protocol: Under the **Graphical** or **Tabular View** the default
- 687           profile should be 2 holding stages and a 2-step cycling protocol. Each stage will have 3 user-
- 688           editable text boxes. The first box value represents the Ramp Rate (%) for that stage, the
- 689           second box value represents the temperature (°C) and the third box value represents the
- 690           time (minutes:seconds).
- 691
- 692                   i. Make the following changes to the default Thermocycler protocol:
- 693                   1. Stage 1 First **Holding Stage**
- 694                           a. Ramp Rate: 100%
- 695                           b. Temp: 55
- 696                           c. Time: 5:00
- 697                   2. Step 1 Second **Holding Stage**.
- 698                           a. Ramp Rate: 100%
- 699                           b. Temp: 60
- 700                           c. Time: 5:00
- 701                   3. Highlight the second **Holding Stage** and select the **Add Stage** button. In the drop
- 702                   down menu select **Holding**
- 703                   4. Step 1 **Third Holding Stage**
- 704                           a. Ramp Rate: 100%
- 705                           b. Temp: 65
- 706                           c. Time: 5:00
- 707                   5. First **2-Step Cycling Stage**
- 708                           a. Number of cycles: 10
- 709                           b. Do NOT check Enable Auto Delta
- 710                           c. Step 1
- 711                                   i. Ramp Rate: 100%
- 712                                   ii. Temp: 92
- 713                                   iii. Time: 0:05
- 714                           d. Step 2
- 715                                   i. Ramp Rate: 100%
- 716                                   ii. Temp: 57
- 717                                   iii. Time: 0:40
- 718                                   iv. Turn data collection "Off" by selecting the **Data Selection** button
- 719                                   at the bottom of the step.
- 720                   6. Highlight step 2 and select the **Add Stage** button. In the drop down menu select
- 721                   **Cycling**
- 722                   7. Second 2-Step **Cycling Stage**
- 723                           a. Number of cycles: 30
- 724                           b. Do NOT check Enable Auto Delta
- 725                           c. Step 1

- 726 i. Ramp Rate: 100%
- 727 ii. Temp: 92
- 728 iii. Time: 0:05
- 729 d. Step 2
- 730 i. Ramp Rate: 100%
- 731 ii. Temp: 57
- 732 iii. Time: 0:40
- 733 iv. Ensure the data collection has been turned “On” for this step
- 734 (default setting)
- 735 8. If a wrong stage is added the stage can be removed by pressing the **Undo “Add**
- 736 **Stage”** button immediately after adding the stage or highlight the stage between
- 737 the vertical lines and select the **Delete Selected** button
- 738 5. Set threshold for each analyte
- 739 a. Select the **Analysis** tab in the upper left menu.
- 740 b. Select **Analysis Settings** button in the top right corner.
- 741 c. Highlight SARS-CoV-2 and deselect the **Use Default Settings** box. De-select **Automatic**
- 742 **Threshold** and change threshold to 75,000. Leave **Automatic Baseline** selected.
- 743 d. Highlight PRC and de-select the **Use Default Settings** box. De-select **Automatic Threshold**
- 744 and change threshold to 10,000. Leave **Automatic Baseline** selected.
- 745 e. At the bottom of the box select **Apply Analysis Settings** button
- 746

Target	Threshold	Baseline Start	Baseline End
SARS-CoV-2	75,000	Auto	Auto
PRC	10,000	Auto	Auto

- 747
- 748 i. Save the new protocol as a template for future use.
- 749 i. At the top of the screen select the drop down menu next to **Save**
- 750 ii. Choose **Save as Template**
- 751 iii. Save in an appropriate folder
- 752 iv. **File name:** ‘Lyra SARS-CoV-2’
- 753 v. **Save as type:** ‘Experiment Document Template files (\*.edf)’
- 754 vi. Exit the software.
- 755

**Applied Biosystems 7500 Standard Thermocycler Test Procedure**

- 756
- 757 1. Launch the 7500 Standard software v2.06 package.
- 758 2. The **Quick Startup document** dialog window will open.
- 759 3. Click on **Create a new document.**
- 760 4. Most of the following should be the default setting. If not, change accordingly.

<b>Assay:</b>	Standard Curve (Absolute Quantitation)
<b>Container:</b>	96-Well Clear
<b>Template:</b>	Lyra SARS-CoV-2
<b>Run Mode:</b>	Fast 7500
<b>Operator:</b>	<i>your operator name</i>
<b>Comments:</b>	SDS v1.4
<b>Plate Name:</b>	<b>YYMMDD-</b> Lyra SARS-CoV-2

- 761 5. Set Up Sample Plate
- 762 a. Under the **Setup** and **Plate** tabs the plate setup will appear.

- 763 b. Select all wells that will contain sample, right-click and select the **Well Inspector** from the drop-  
764 down menu. When the **Well Inspector** pop-up window opens, select the detectors for SARS-CoV-  
765 2 and PRC.
- 766 c. Use the **Well Inspector** to enter the sample names. Patient IDs can be entered in the Well Inspector  
767 window. However, it is recommended that this is done prior to re-suspending the lyophilized  
768 master mix, post run or using the import function to minimize the time the PCR reactions will sit  
769 at room temperature prior to starting the run.
- 770 d. Save the run as **YYMMDD- Lyra SARS-CoV-2.sds**.
- 771 e. A window will open asking for the “Reason for change of entry”. Enter “**Setup**” and any other  
772 comments relevant to the run.
- 773 6. Starting the PCR
- 774 a. Select the **Instrument** tab.
- 775 b. Insert the 96 well PCR plate into the machine.
- 776 c. Under **Instrument Control**, select the **Start** button to initiate the run.
- 777 7. Post PCR
- 778 **IMPORTANT:** When the run is finished press OK.
- 779 a. Analyze the data by pressing the “**Analyze**” button in the top menu and save the file.
- 780 b. Save the file by pressing **Save Document** in the task bar. A window will open asking for the  
781 “Reason for change of entry”.
- 782 c. Enter “**Data analysis post run**” and any other comments relevant to the run.
- 783

## 784 Bio-Rad CFX96 Touch Thermocycler Programming Procedure

785 Refer to User Manual Part Number 10010424 Rev D for additional information.

786

### 787 Programming Instructions:

- 788 1. Launch the CFX96 Touch software package
- 789 2. In the **Startup Wizard** pop-up window **Select instrument** to be **CFX96** from the drop down menu
- 790 3. Under **Select Run Type** press the **User-defined** button
- 791 4. Create a new thermocycler protocol by selecting **Create New** from the **Run Setup** window
- 792 5. Make the following changes to the cycling conditions in the **Protocol Editor**:
- 793 a. Change the **Sample Volume** to **20ul**
- 794 b. Under **Tools** in the top left toolbar select **Run Time Calculator** and check **96 Wells-All Channels**
- 795 c. **Step 1** (Hold)
- 796 i. Reps: 1
- 797 ii. Temp: 55C
- 798 iii. Time: 5:00
- 799 d. **Step 2** (Hold)
- 800 i. Reps: 1
- 801 ii. Temp: 60C
- 802 iii. Time: 5:00
- 803 e. **Step 3** (Hold)
- 804 i. Reps: 1
- 805 ii. Temp: 65C
- 806 iii. Time: 5:00
- 807 iv. Remove the plate read from this stage by selecting the **Remove Plate Read** button on the  
808 lower left
- 809 f. **Step 4** (2-Step Amplification Stage)
- 810 i. Highlight **step 3** and go to the lower left of the window and select **Insert Step** for a total  
811 of 2 times until step 5 is reached (ensure in the upper left of the window the drop-down  
812 menu for **Insert Step** has **After** selected).
- 813 ii. Highlight **step 4** and set as follows:
- 814 1. Temp: 92C

- 
- 815                                   2. Time: 0:05
- 816                                   iii. Highlight **step 5** and set as follows:
- 817                                   1. Temp: 57C
- 818                                   2. Time: 0:40
- 819                                   3. Go to the left of the screen and select **Remove Plate Read** button
- 820                                   iv. Select **step 6**, the **GOTO step**, and change to state **GOTO step 4** and change the times to
- 821                                   repeat to **9**
- 822                                   g. **Step 7** (2-Step Amplification Stage)
- 823                                   i. With step 6 highlighted select **Insert Step** button, on the lower left of the window, for a
- 824                                   total of 2 times (until step 8 is reached)
- 825                                   ii. Highlight **step 7** and set as follows:
- 826                                   1. Temp: 92C
- 827                                   2. Time: 0:05
- 828                                   iii. Highlight **step 8** and set as follows:
- 829                                   1. Temp: 57C
- 830                                   2. Time: 0:40
- 831                                   3. In the left of the window select **Add Plate Read to Step** button
- 832                                   4. Highlight **step 8** and select **Insert GOTO** button on the lower left of the window
- 833                                   iv. Select **step 9**, the **GOTO step**, and change to **GOTO step 7** and times to repeat to **29**
- 834                                   h. Save the new cycling conditions as protocol for future use
- 835                                   i. At the upper left of the screen select the **Save** button
- 836                                   ii. Save in the **ExpressLoad** folder
- 837                                   iii. **Name** the file 'Lyra SARS-CoV-2'
- 838                                   iv. **Save as type** 'Protocol File (\*.prcl)'
- 839                                   v. Select **Save**
- 840                                   vi. Click **Ok** in the protocol editor window
- 841                                   6. Define the plate setup
- 842                                   a. In the **Run Setup** window select the **Plate** tab
- 843                                   b. Under **Express Load** in the drop-down menu select **Quick Plate 96 wells All Channels.pltd**
- 844                                   c. Select the **Edit Selected** button to customize the plate setup
- 845                                   d. In the upper toolbar select **Settings**. The default settings need to be set.
- 846                                   i. **Plate Size** select **96 Wells**
- 847                                   ii. **Plate Type** select **BR Clear**
- 848                                   iii. **Number Convention** select **Scientific Notation**
- 849                                   iv. **Units** select **Copy Number**
- 850                                   e. Leave the **Scan Mode** set to **All Channels** at the top of the window
- 851                                   f. Select the **Select Fluorophores** button on the upper right of the Plate Editor window
- 852                                   i. De-select all default fluorophores
- 853                                   ii. Select **FAM**, and **Cy5** and click Ok
- 854                                   g. In the **Plate Editor** window highlight the whole plate and click the check box in front of all
- 855                                   fluorophores: **FAM** and **Cy5**
- 856                                   h. Select the **Experiment Settings** button in order to define the Targets
- 857                                   i. In the lower left of the **Experiment Settings** window in the **New** box type in **SARS-CoV-2**
- 858                                   and select **Add**
- 859                                   ii. Repeat this for the **PRC**
- 860                                   iii. Select **Ok**
- 861                                   i. In the **Plate Editor** window next to **FAM** in the drop-down menu under **Target Name** select **SARS-**
- 862                                   **CoV-2** and for Cy5 select **PRC**
- 863                                   j. Save the new plate setup for future use
- 864                                   i. At the upper left of the screen select the **Save** button
- 865                                   ii. Save in the **ExpressLoad** folder
- 866                                   iii. **Name** the file 'Lyra SARS-CoV-2 plate'

- 867 iv. **Save as type** 'Plate File (\*.pltd)'
- 868 v. Select **Save**
- 869 vi. Click **Ok** in the **Plate Editor** window
- 870 k. Exit the software
- 871

## 872 **Bio-Rad CFX96 Touch Thermocycler Test Procedure**

### 873 **Analysis Instructions:**

- 874 1. Open the run file that needs to be analyzed
- 875 2. In the upper left select the **Quantification Tab**
- 876 3. On the Amplification curve check the box in front of **Log Scale**
- 877 4. Select **Settings** in the toolbar in the upper left of the screen
  - 878 a. For the **Cq Determination Mode** select **Single Threshold**
  - 879 b. Under the **Baseline Setting** choose **Baseline Subtracted Curve Fit**
  - 880 c. For **Analysis Mode** select **Target**
  - 881 d. Under **Cycles to Analyze** choose 1-30 and then click **Ok**
  - 882 e. The baseline cycles and the threshold for each target need to be set
    - 883 i. Ensure that only the **SARS-CoV-2 box** is checked in the amplification plot
    - 884 ii. Go up to **Settings** in the toolbar and select **Baseline Threshold**
      - 885 1. At the top of the box select **Auto Calculated** for the **Baseline Cycles**
      - 886 2. For the **Single Threshold** at the bottom of the box select **User Defined**
        - 887 a. Set this to **164**
        - 888 b. Select **Ok**
    - 889 iii. **Uncheck** the **SARS-CoV-2 box** and **check** the **PRC box** in the amplification plot
    - 890 iv. Go up to **Settings** in the toolbar and select **Baseline Threshold**
      - 891 1. At the top of the box select **Auto Calculated** for the **Baseline Cycles**
      - 892 2. For the **Single Threshold** at the bottom of the box select **User Defined**
        - 893 a. Set this to **100**
        - 894 b. Select **Ok**
  - 895 5. Exit the software
  - 896

## 897 **Qiagen Rotor-Gene Q Programming Instructions**

898 Refer to User Manual Part Number 1065453EN for additional information.

### 899 **Programming Instructions:**

- 900 1. Launch the Rotor-Gene Q software package
- 901 2. In the **New Run** pop-up window select the **Advanced** tab on the top of the screen
- 902 3. Select **Empty Run** and then **New** on the lower right of the pop-up window to start the **Advanced Run Wizard**
  - 903 a. Select the appropriate rotor size in the **Advanced Run Wizard** on the upper left of the screen
  - 904 b. Check the box that states the **Locking Ring is Attached** and select **Next**
  - 905 c. Leave the **Operator** and **Notes** sections empty
  - 906 d. Enter **20ul** as the **Reaction Volume** in the lower left of the screen
  - 907 e. For the **Sample Layout** choose **1, 2, 3...** and then select **Next**
  - 908 f. Under **Channel Setup** select **Create New** to enter information for each detector
    - 909 i. Under **Name** enter **SARS-CoV-2**
    - 910 ii. **Source** select 470nm
    - 911 iii. **Detector** select 510nm
    - 912 iv. Do not adjust the default **Gain** setting of 7 as this will be set in a later step
    - 913 v. Select **OK**
  - 914 g. Repeat the step above by selecting **Create New**
    - 915 i. Under **Name** enter **PRC**
    - 916 ii. **Source** select 625nm
    - 917 iii. **Detector** select 660nm
    - 918 iv. Do not adjust the default **Gain** setting of 7 as this will be set in a later step

- 
- 919 v. Select **OK**
- 920 h. Select the **Edit Profile** button in the middle of the window to setup a cycling profile
- 921 i. In the **Edit Profile** window go to the upper left of the screen to **New** and in the drop-down
- 922 menu select **Cycling**. A hold and three step cycling stage should appear.
- 923 ii. Modify the hold stage to have a temperature at **55°C** and a time of **5:00 minutes**
- 924 iii. Select the **Insert After** button in the middle of the pop-up window and then select **New**
- 925 **Hold at Temperature**
- 926 iv. Modify the second hold stage to have a temperature at **60°C** and a time of **5:00 minutes**
- 927 v. Select the **Insert After** button in the middle of the pop-up window and then select **New**
- 928 **Hold at Temperature** to insert a third hold stage
- 929 vi. Modify the third hold stage to have a temperature at **65°C** and a time of **5:00 minutes**
- 930 vii. Highlight the first **cycling stage** and modify it as follows:
- 931 1. This cycle repeats **10** time(s)
- 932 2. Select **Timed Step** from the drop-down menu in the middle left of the screen
- 933 3. Do not select **Long Range** or **Touchdown** on the left of the screen
- 934 4. The first step:
- 935 a. **92°C**
- 936 b. **5 seconds**
- 937 c. **Not Acquiring**
- 938 5. Select step two and set as follows:
- 939 a. **57°C**
- 940 b. **40 seconds**
- 941 c. **Not Acquiring**
- 942 6. Highlight step three and delete it by selecting the “-” button in the middle of the
- 943 window
- 944 7. Select the **Insert After** button in the middle of the pop-up window and then
- 945 select **New Cycling**
- 946 viii. Highlight the second **cycling stage** and modify it as follows:
- 947 1. This cycle repeats **30** time(s)
- 948 2. Select **Timed Step** from the drop-down menu in the middle left of the screen
- 949 3. Do not select **Long Range** or **Touchdown** on the left of the screen
- 950 4. The first step:
- 951 a. **92°C**
- 952 b. **5 seconds**
- 953 c. **Not Acquiring**
- 954 5. Select step two and set as follows:
- 955 a. **57°C**
- 956 b. **40 seconds**
- 957 c. Select **Acquiring to Cycling A**
- 958 i. **Under Acquiring Channels** highlight the default channel name
- 959 (Green) and select the < button to move it over to the
- 960 **Available Channels** list
- 961 ii. In the **Available Channels** list select **SARS-CoV-2** and select the
- 962 > button to move it over to the **Acquiring Channels** list
- 963 iii. Repeat the step above for the **PRC** and then select **OK**
- 964 6. Highlight step three and delete it by selecting the “-” button in the middle of the
- 965 window
- 966 ix. In the **Edit Profile** window select **OK**
- 967 i. In the **New Run Wizard** window select **Gain Optimisation**
- 968 i. In the middle of the **Auto-Gain Optimisation Setup** window select the drop-down menu
- 969 under **Channel Settings** and select **SARS-CoV-2**.
- 970 ii. Select the **Add** button on the right

- 971 1. In the **Auto-Gain Optimisation Channel Settings** window ensure that the SARS-
- 972 CoV-2 **Tube Position** is set to **1**. This requires that a positive control, containing
- 973 SARS-CoV-2 and PRC, be tested with each PCR run and placed in the first tube.
- 974 Failure to do so may cause the gain to be incorrectly set.
- 975 2. Leave the **Target Sample Range** and the **Acceptable Gain Range** set to the
- 976 defaults, 5-10FI and -10 to 10 respectively.
- 977 3. Select **OK**
- 978 4. Repeat steps 3. j. ii. 1-3. for the **PRC**
- 979 **iii.** In the **Auto-Gain Optimisation Setup** window check the box next to **Perform**
- 980 **Optimisation Before 1<sup>st</sup> Acquisition**
- 981 **iv.** Select **Close**
- 982 **j.** In the **New Run Wizard** window select the **Next** button
- 983 **k.** Save the new protocol as a template for future use
- 984 **i.** On the bottom right of the window select the **Save Template** button
- 985 **ii.** **Save In:** C:\Program Files\Rotor-Gene Q Software\Templates
- 986 **iii.** **File name:** 'Lyra SARS-CoV-2'
- 987 **iv.** **Save as type:** 'Template (\*.ret)'
- 988 **l.** Exit the software
- 989

## 990 Qiagen Rotor-Gene Q Test Run

### 991 Analysis Instructions:

- 992 1. In the New Run Wizard load the SARS-CoV-2 Template.
- 993 2. Press Start.
- 994 3. Open the run file that needs to be analyzed
- 995 4. In the upper menu toolbar select the **Analysis** button
- 996 **a.** Select **Quantitation**, then **Cycling A. SARS-CoV-2**, and **Show**
- 997 **b.** The threshold needs to be set for SARS-CoV-2
- 998 **i.** In the far right bottom of the screen under **CT Calculation** enter **0.03** for the **SARS-CoV-2**
- 999 **Threshold**
- 1000 **ii.** In the **Eliminate Cycles before** box ensure the default of **1** is entered
- 1001 **iii.** Ensure the amplification graph is set to **Log Scale** (toggle button on the bottom left of the
- 1002 graph states Linear Scale or Log Scale)
- 1003 **c.** Select **Quantitation**, then **Cycling A. PRC**, and **Show**
- 1004 **d.** The threshold needs to be set for PRC
- 1005 **i.** In the far right bottom of the screen under **CT Calculation** enter **0.05** for the **PRC**
- 1006 **Threshold**
- 1007 **ii.** In the **Eliminate Cycles before** box ensure the default of **1** is entered
- 1008 **iii.** Ensure the amplification graph is set to **Log Scale** (toggle button on the bottom left of the
- 1009 graph states Linear Scale or Log Scale)
- 1010

## 1011 Roche LightCycler 480 Instrument II Programming Instructions

1012 Refer to User Manual Part Number 05152062001 0208 for additional information.

### 1014 Creating a LightCycler 480 II Assay Run Template

- 1015 1. Launch the LightCycler (LC) 480 software package
- 1016 2. The **Detection Format** must be established to specify the channels in which fluorescence will be read
- 1017 **a.** Select **Tools** in the startup screen in the lower right of the screen
- 1018 **b.** Select **Detection Formats** then choose **New**
- 1019 **c.** Name the format Lyra® SARS-CoV-2
- 1020 **d.** In the **Filter Combination Selection** window select 465-510 and 618-660

- 
- 1021 e. In the **Selected Filter Combination List** window under name type in SARS-CoV-2 for 465-510 and  
1022 PRC for 618-660
- 1023 f. Leave all default setting values to 1 under Melt Factor, Quant Factor, and Max Integration Time
- 1024 g. Select **Close** to save the new detection format and return to startup screen
- 1025 h. To access this newly created **Detection Format**, the LightCycler 480 software must be closed, then  
1026 reloaded
- 1027 3. After closing and reloading the software select **White Plates** and **New Experiment** under Experiment  
1028 Creation window
- 1029 4. On the next screen select “Lyra® SARS-CoV-2” from the pull-down menu under **Detection Formats**
- 1030 5. Enter **20ul** as the **Reaction Volume** in the upper right of the screen
- 1031 6. Enter the names for each of the RT-PCR programs
- 1032 a. Under **Program Name** enter **Stage 1**, under **Cycles** enter **1**, and in **Analysis Mode** select **none**
- 1033 b. Select the “+” icon to add a program
- 1034 c. Name the next program **Stage 2**, under **Cycles** enter **1**, and in the **Analysis Mode** select **none**
- 1035 d. Select the “+” icon to add a program
- 1036 e. Name the next program **Stage 3**, under **Cycles** enter **1**, and in the **Analysis Mode** select **none**
- 1037 f. Select the “+” icon to add a program
- 1038 g. Name the next program **Stage 4**, under **Cycles** enter **40**, and in the **Analysis Mode** select  
1039 **quantification**
- 1040 7. Set the RT-PCR cycling times and temperatures
- 1041 a. Highlight **Stage 1** under **Program Name** and change **Stage 1 Temperature Targets** as follows:
- 1042 i. **Target (°C)** set to **55**
- 1043 ii. **Acquisition Mode** select **none**
- 1044 iii. **Hold (hh:mm:ss)** set to **5:00**
- 1045 iv. **Ramp Rate (°C/s)** to 4.4
- 1046 v. **Sec Target (°C)**, **Step Size (°C)**, and **Step Delay (cycles)** will be left at 0 for stages 1-4.
- 1047 b. Highlight **Stage 2** under **Program Name** and change **Stage 2 Temperature Targets** as follows:
- 1048 i. **Target (°C)** set to **60**
- 1049 ii. **Acquisition Mode** select **none**
- 1050 iii. **Hold (hh:mm:ss)** set to **5:00**
- 1051 iv. **Ramp Rate (°C/s)** to 4.4
- 1052 c. Highlight **Stage 3** under **Program Name** and change **Stage 3 Temperature Targets** as follows:
- 1053 i. **Target (°C)** set to **65**
- 1054 ii. **Acquisition Mode** select **none**
- 1055 iii. **Hold (hh:mm:ss)** set to **5:00**
- 1056 iv. **Ramp Rate (°C/s)** to 4.4
- 1057 d. Highlight **Stage 4** under **Program Name** and change **Stage 4 Temperature Targets** as follows:
- 1058 i. The first step:
- 1059 1. **Target (°C)** set to **92**
- 1060 2. **Acquisition Mode** select **none**
- 1061 3. **Hold (hh:mm:ss)** set to **0:05**
- 1062 4. **Ramp Rate (°C/s)** to 4.4
- 1063 ii. Select the “+” icon to add a step and set the second step:
- 1064 1. **Target (°C)** set to **57**
- 1065 2. **Acquisition Mode** select **single**
- 1066 3. **Hold (hh:mm:ss)** set to **0:40**
- 1067 4. **Ramp Rate (°C/s)** to 2.2
- 1068 8. Save the new protocol as a run template for future use.
- 1069 a. In the lower left corner of the screen select the pull-down menu next to the **Apply Template** button
- 1070 b. Choose **Save As Template**
- 1071 c. Select the **Templates Folder**
- 1072 d. Highlight **Run Templates Folder**



- 1073 e. Name the template Lyra® SARS-CoV-2 run template and click the “check” button  
1074 9. Exit the software.  
1075

### 1076 **Creating a Roche LightCycler 480 II Assay Test Procedure**

- 1077 1. Load the Lyra SARS-CoV-2 run template.  
1078 2. Press Start.  
1079 3. The analysis template can only be established after the initial experiment has completed  
1080 4. On the Lyra® SARS-CoV-2 run select the **Analysis** button in the module bar  
1081 a. Choose **Abs Quant/Fit Points**  
1082 b. In the **Create New Analysis** pop-up window select your pre-defined subset from the **subset** drop  
1083 down menu and then select the “check” button  
1084 c. Set the **Background** to 2-10 for all analytes  
1085 i. Set **Min Offset** to 1  
1086 ii. Set **Max Offset** to 9  
1087 d. In the center bottom of the screen ensure that **Color Compensation** is off for all analytes  
1088 e. Leave the default settings as **First Cycle 1** and **Last Cycle 40**  
1089 5. At the top middle of the screen select **Noise Band**  
1090 6. Choose the pull-down menu next to the **Noise Band** button and select **Noise Band Fluorescence**  
1091 7. For each analyte under the **Filter Comb** button, set the noise band as follows:  
1092 a. SARS-CoV-2 set to 1.95  
1093 b. PRC set to 1.4619  
1094 8. Choose **Calculate** in the bottom left of the screen  
1095 9. Save the new analysis protocol as a template for future use  
1096 a. In the lower left corner of the screen select the pull-down menu next to the **Apply Template** button  
1097 b. Choose **Save As Template**  
1098 c. Select the **Templates Folder**  
1099 d. Highlight **Analysis Templates Folder**  
1100 e. Name the template Lyra® SARS-CoV-2 analysis template and click the “check” button  
1101 10. Create a report  
1102 a. Select the **Save** icon on the global action bar on the right side of the screen  
1103 b. Choose the **Report** button on the module bar on the left of the screen  
1104 c. Select the appropriate settings and press the **Generate** button  
1105 11. To apply an Analysis Template to subsequent runs  
1106 a. Once the run has finished select the **Analysis** button in the module bar  
1107 b. Choose **Abs Quant/Fit Points**  
1108 c. In the **Create New Analysis** pop-up window select your pre-defined subset from the **subset** drop  
1109 down menu and then select the “check” button  
1110 d. Select the **Apply Template** button on the far left of the screen and choose the Lyra® SARS-CoV-2  
1111 analysis template from the **Analysis Templates Folder**  
1112 e. Select yes in the pop-up window  
1113 12. Interpretation of results (See Table 4)  
1114

### 1115 **Thermo Fisher QuantStudio 7 Pro Programming Instructions**

1116 Refer to User Manual Part Number 4489822 Revision A for additional information.

1117

### 1118 **Thermo Fisher QuantStudio 7 Pro Test Run Programming Instructions:**

- 1119 1. Open the Design and Analysis Software  
1120 2. Select the “SET UP PLATE” option  
1121 3. From the side bar on the screen, select the following properties to filter:  
1122 a. Instrument – QuantStudio 7 Pro  
1123 b. Block – 96-Well 0.2 mL

- 
- 1124 c. Run Mode – Fast
- 1125 d. Analysis options are left blank
- 1126 4. From the plate selections present on the screen, select the System Template “PCR Only” and
- 1127 the system will automatically navigate to the “Run Method” tab
- 1128 5. Run Method
- 1129 a. Change the Reaction Volume to 20.0 uL
- 1130 b. The temperature of the enabled heated cover will remain at 105.0 degrees C
- 1131 c. Scroll over the Hold stage present in the cycling parameters and
- 1132 addition/subtraction buttons will become visible at both the top and bottom of the
- 1133 first stage.
- 1134 d. Left click the right addition button at the top and a list of Stage choices will become
- 1135 visible. Scroll down and choose Hold.
- 1136 e. Repeat the previous steps so there are three Hold stages present in the cycling
- 1137 parameters.
- 1138 f. Scroll over to the PCR stage and addition/subtraction buttons will become visible at
- 1139 both the top and bottom. Left click the right addition button at the top and a list of
- 1140 Stage choices will become visible. Scroll down and choose PCR.
- 1141
- 1142 g. Going back to the first stage enter the following parameters:
- 1143 i. Stage 1 Hold
- 1144 1. 2.63 ramp rate
- 1145 2. 55°C
- 1146 3. 5 minutes
- 1147 ii. Stage 2 Hold
- 1148 1. 2.63 ramp rate
- 1149 2. 60°C
- 1150 3. 5 minutes
- 1151 iii. Stage 3 Hold
- 1152 1. 2.63 ramp rate
- 1153 2. 65°C
- 1154 3. 5 minutes
- 1155 iv. Stage 4 PCR
- 1156 1. Step 1:
- 1157 a. 2.63 ramp rate
- 1158 b. 92°C
- 1159 c. 5 seconds
- 1160 2. Step 2:
- 1161 a. 2.32 ramp rate
- 1162 b. 57°C
- 1163 c. 40 seconds
- 1164 d. Click on the camera icon under Step 2. A window will pop
- 1165 up asking for confirmation to turn off data collection during
- 1166 this step. Click “Ok”.
- 1167 v. Located at the bottom of Stage 4 PCR change the number of cycles to 10
- 1168 vi. Stage 5 PCR
- 1169 1. Step 1:
- 1170 a. 2.63 ramp rate

- 1171 b. 92°C  
1172 c. 5 seconds  
1173 2. Step 2:  
1174 a. 2.32 ramp rate  
1175 b. 57°C  
1176 c. 40 seconds  
1177 d. Ensure the camera icon image is bold/on for data collection  
1178 during the 30 cycles of Stage 5, Step 2.  
1179 vii. Located at the bottom of Stage 4 PCR change the number of cycles to 30  
1180 h. Scroll up and choose the “Plate Setup” tab near the top of the screen.  
1181 6. Plate Setup  
1182 a. Change the Passive Reference to “NONE”  
1183 b. On the lower right side of screen, ensure the Targets Tab is chosen then highlight  
1184 and press the addition button to add “Target 1”. Press again to add “Target 2”  
1185 c. Click on the “Target 1” box and change the name to CoV-2.  
1186 d. Click the associated reporter box below the Reporter tab and, from the drop down  
1187 menu, choose FAM.  
1188 e. Click on the “Target 2” box and change the name to PRC.  
1189 f. Click the associated reporter box below the Reporter tab and, from the drop down  
1190 menu, choose CY5.  
1191 g. Highlight the “Actions” button located in the upper right side of the screen and  
1192 press the drop down button. In the drop down menu choose “Analysis Setting”  
1193 h. Under Analysis Setting, disable the following for all targets:  
1194 i. Use Default Column  
1195 ii. Auto Threshold Column  
1196 iii. Auto Baseline Column  
1197 iv. The Baseline Start and Baseline End should default to 3 and 15  
1198 i. Under “Threshold” click on the box associate with the CoV target and enter 70000.  
1199 j. Under “Threshold” click on the box associated with the PRC target and enter 20000  
1200 k. Click “Save”  
1201 l. Navigate back to the “Actions” button and press the drop down button, choosing  
1202 “Save As”. This will save your template to a location of choice. Save the template as  
1203 “Lyra SARS Cov-2 Assay”.  
1204

## 1205 **Creating a Thermo Fisher QuantStudio 7 Pro Test Procedure**

1206

1207 **Note:** These instructions are based upon the user not having the QuantStudio 7 Pro Real-Time PCR  
1208 instrument and the Applied Biosystems Design and Analysis 2.2 software connected. The user must  
1209 open the Lyra SARS CoV-2 template created previously with the software and save any newly  
1210 created sample run template onto a USB and transfer the template to the instrument.  
1211 For connectivity related to the software and the instrument please contact your Thermo Fisher  
1212 QuantStudio representative.  
1213

- 1214 1) Open the Lyra SARS CoV-2 Assay Template previously generated.  
1215 2) Click on the Plate Setup Tab located near the top of the screen.  
1216 3) On the right side of the screen ensure the “Samples” tab is highlighted and press the  
1217 addition button to add the number of samples being tested.

- 
- 1218 4) Click on the “Sample 1” box to rename the sample. Repeat this step for all subsequent  
1219 samples being entered.
- 1220 5) Click the well located in the plate map then check the box next to the sample name from the  
1221 right side bar to associate the name to the well.
- 1222 a. User also has the option to highlight the well location in the plate map and click on  
1223 the “Enter sample” box. Enter the sample ID and press tab to continue to the next  
1224 well in the plate map. This will automatically load the sample name into the  
1225 sidebar.
- 1226 6) Once samples names have been entered, the wells may be highlighted by left clicking the  
1227 mouse over starting well and dragging the mouse across all wells associated in run. The  
1228 targets are then chosen by clicking the check boxes next to each target in the side bar.
- 1229 7) Click on the Actions button located top right of the screen and choose “Save As” in the  
1230 dropdown menu.
- 1231 a. A pop-up window will appear directing the user to title the file according to  
1232 information pertaining to the sample run and the location of the file to be saved.
- 1233 b. Save the newly named (.edt) run file to a USB that is inserted into the computer.
- 1234 8) Transfer the USB to the port on the front of the instrument.
- 1235 9) From the options on the instrument’s screen press “Load plate file”. The QuantStudio 7 Pro  
1236 is a touchscreen device.
- 1237 10) From the “Run Queue” screen, press “USB drive” on the right side. This will bring up any  
1238 plate files saved on the USB.
- 1239 11) Press the plate file associated with the run to be performed.
- 1240 12) A new window will appear requesting location of results once the run is complete.
- 1241 a. Press the “USB drive Connected” if the icon is not already highlighted and press  
1242 “Done”.
- 1243 13) Centrifuge the 96-well sample plate to ensure all liquid is toward the bottom of each well.
- 1244 a. Ensure the centrifuge is properly balanced.
- 1245 b. Gently pull the plate from the centrifuge to ensure all liquids remain at the bottom  
1246 of the wells.
- 1247 14) Press the double-arrow icon located at the top right sided corner of the screen on the  
1248 instrument.
- 1249 a. The instrument drawer will open from the front.
- 1250 15) Place the centrifuged plate into the plate holder ensuring proper orientation of the plate.
- 1251 a. A1 well should be in the position of the top left corner
- 1252 b. The plate will appear slightly suspended above the block due to two silicone strips  
1253 above and below this plate. This is to be expected and the instrument lid will press  
1254 the plate down once the drawer has closed.
- 1255 16) Press “Start Run” on the screen of the instrument.
- 1256 a. A pop-up window will appear asking the user to confirm the plate has been loaded.
- 1257 b. If the plate has been loaded, press “Start Run” again or press “Open Drawer” to  
1258 place the plate into the block and then press “Start Run”  
1259

1260



M120 – Lyra SARS-CoV-2 Assay kit



**Quidel Corporation**

2005 East State Street, Suite 100

Athens, OH 45701 USA

**quidel.com**

1261

1262

---

**GLOSSARY**



Intended use



Catalog number



Contents / Contains



Contains sufficient for 96 determinations



Control



Batch code



Use by



Consult e-labeling instructions for use



Manufacturer



Temperature limitation