



# Solana<sup>®</sup> Influenza A+B Control Set

## Positive and Negative

QUIDEL



### INTENDED USE

The Solana Influenza A+B Control Set is intended to be used as full process controls in the Solana Influenza A+B Assay. This unassayed control set for *in vitro* diagnostic use contains a buffered solution of synthetic protected RNA sequences for influenza A and influenza B primer and probe targets. The negative control is an influenza A and B free matrix.

### INSTRUCTIONS FOR USE

The Solana Influenza A+B controls should be used at the same volume as unknown samples with the Solana Influenza A+B Assay.

### MATERIALS PROVIDED

Cat. #M122

Control Kit (2 vials) – Store at 2°C to 8°C

#	Component	Quantity
<b>CONTROL</b> +	<b>Influenza A+B Positive Control</b> Part M5126	1 vial/kit 1.0 mL per vial
<b>CONTROL</b> -	<b>Negative Control</b> Part M5203	1 vial/kit 1.0 mL

### WARNINGS AND PRECAUTIONS

- Treat all specimens, samples and controls as potentially infectious. Follow universal precautions when handling samples, this kit, and its contents.
- Avoid contact with skin and eyes.
- Use micropipettes with an aerosol barrier or positive displacement tips for all procedures.
- Avoid microbial and cross contamination of the kit reagents. Follow Good Laboratory Procedures.
- Controls contain 0.09% Sodium Azide. Do not empty these controls into drains. Dispose of the containers in a safe way.
- Proper sample collection, storage, and transport are essential for correct results.<sup>1</sup>
- Store assay reagents as indicated on their individual labels.
- Testing should be performed in an area with adequate ventilation.
- Dispose of containers and unused contents in accordance with Federal, State and Local regulatory requirements.
- Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this kit.
- Wash hands thoroughly after handling.
- For additional information on hazard symbols, safety, handling and disposal of the components within this kit, please refer to the Safety Data Sheet (SDS) located at [quidel.com](http://quidel.com).

### STORAGE AND HANDLING OF KIT REAGENTS

Store the kit at 2°C to 8°C until the expiration date listed on the labels.

## QUALITY CONTROL

The Solana Influenza A+B Control Set can be used for quality control testing. Each laboratory should establish their own frequency of QC testing based on applicable local laws, regulations and standard good laboratory practice.

## LIMITATIONS

- For *in vitro* diagnostic use.
- This control has been studied using a molecular testing methodology; it is not intended for use with other methodologies.

## CUSTOMER AND TECHNICAL SUPPORT

To place an order or for technical support, please contact a Quidel Representative at 800.874.1517 (in the U.S.) or 858.552.1100 (outside the U.S.), Monday through Friday, from 8:00 a.m. to 5:00 p.m., Eastern Time. Orders may also be placed by fax at 740.592.9820. For e-mail support, contact [customerservice@quidel.com](mailto:customerservice@quidel.com) or [technicalsupport@quidel.com](mailto:technicalsupport@quidel.com).

For services outside the U.S., please contact your local distributor. Additional information about Quidel, our products, and our distributors can be found on our website [quidel.com](http://quidel.com).

## REFERENCES

1. CLSI MM13-A: Guidance for the Collection, Transport, Preparation and Storage of Specimens for Molecular Methods 25:31, 2005.

**REF** M122 – Solana Influenza A+B Control Set

**IVD**



**EC REP**

MDSS GmbH  
Schiffgraben 41  
30175 Hannover,  
Germany



**Quidel Corporation**  
2005 East State Street, Suite 100  
Athens, OH 45701 USA  
[quidel.com](http://quidel.com)

PIM122000EN00 (09/17)

## GLOSSARY

---

**REF**

Catalogue number



CE mark of conformity

---

**EC REP**

Authorized representative  
in the European Community

**LOT**

Batch code

---



Use by



Manufacturer

---



Temperature limitation



Intended use

---



Biological risks



Consult e-labeling  
instructions for use

---

**CONT**

Contents/Contains

**IVD**

For *in vitro* diagnostic use

---

**CONTROL +**

Control positive

**CONTROL -**

Control negative

---