



QUIDEL®

Technical Bulletin

SARS-CoV-2 Variant Detection

A study was performed to determine detection of various SARS-CoV-2 recombinant mutations across the on-market Quidel antigen assays.

Each recombinant mutation, as specified by the proposed variants and mutations to the nucleoprotein, was diluted to an initial target concentration and tested on the Sofia® SARS Antigen FIA to verify that it was detectable. Once detected, serial dilutions of each recombinant mutation were subsequently formulated and tested with each of the Quidel antigen assays listed in the table below. Each level was tested with a total number of three (3) replicates. All eight (8) recombinant mutations were detected by each of the Quidel antigen assays.

Variant Strain	Place of Origin	NP Mutation	Sofia SARS Antigen	Sofia 2 Flu + SARS Antigen	QuickVue® SARS Antigen	QuickVue At-Home COVID-19
Wild type	Wuhan	N/A	✓	✓	✓	✓
B.1.1.7	UK	D3L, S235F	✓	✓	✓	✓
B.1.1.7	UK	D3L, R203K, S235F	✓	✓	✓	✓
B.1.1.7	UK	D3L, R230K, G240R, S235F	✓	✓	✓	✓
B.1.1.7	UK	S194L	✓	✓	✓	✓
501Y.V2	South Africa	T205I	✓	✓	✓	✓
B.1.1.7, 501Y.V2	UK/South Africa	D3L/T205I/S235F	✓	✓	✓	✓
P.1 (501Y.V3)	Brazil	R203K, G204R	✓	✓	✓	✓
P.1 (501Y.V3)	Brazil	I292T	✓	✓	✓	✓
P.1 (501Y.V3)	Brazil	A119S-R203K-G204R-M234	✓	✓	✓	✓
I/A3i	India	P13L	✓	✓	✓	✓

If you have any questions regarding the use of this product or to report a product problem, please contact Quidel Technical Support at 1.800.874.1517 (in the U.S.) or technicalsupport@quidel.com. If outside the U.S., further information can be obtained from your distributor, or directly from Quidel at one of the numbers listed below. Reference quidel.com to see more options for Support.

Country	Phone	E-Mail Address
Europe, Middle East and Africa	+353 (91) 412 474 (main) 0 1800 200441 (toll free)	emeatechnicalsupport@quidel.com
Austria	+43 316 231239	
France	0 (805) 371674	
Germany	+49 (0) 7154 1593912	
Netherlands	0 800 0224198	
Switzerland	0 800 554864	
United Kingdom	0 800 3688248	
Italy	+39 (800) 620 549	
North America, Asia-Pacific, Latin America	858.552.1100	technicalsupport@quidel.com
Canada	437.266.1704 (main) 888.415.8764 (toll free)	technicalsupport@quidel.com
China	0400 920 9366 or +86 021 3217 8300	chinatechnicalservice@quidel.com

You may also visit our website at quidel.com for information on Quidel's line of Rapid Diagnostics, Molecular Diagnostics, Cell Culture and Specialty Products (Bone Health and Autoimmune & Complement). Other product information available on our website includes: CPT codes, CLSI procedure guides, SDS, and Package Inserts.

TB10181701EN00 (04/21)