

QuickVue® RSV 10 Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue RSV 10 test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current package insert instructions and/or procedure card included in the kit.

1. How often does the manufacturer recommend that external controls be run on the kit?

- a) Once per kit
- b) Once per each new shipment of kits
- c) Once per day
- d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested

2. At what temperature do the kits need to be stored?

- a) -20°C (frozen)
- b) 15-30°C (room temp.)
- c) 2-30°C (refrigerated or room temp.)
- d) 2-8°C (refrigerated)

3. To guarantee the performance characteristics indicated in the package insert, which swab should be used when collecting a nasopharyngeal swab specimen?

- a) Cotton-tipped swabs
- b) Any type of swabs
- c) Any nasal swab
- d) The swabs provided in the kit

4. What is the read time for this kit?

- a) 10 minutes
- b) 15 minutes
- c) 20 minutes
- d) 8 minutes

5. Can the test be read after the read time?

- a) Yes, any time after the read time
- b) Yes, up to 5 minutes after the read time
- c) No
- d) Yes, up to 10 minutes after the read time

6. How long after adding the Reagent solution to the Reagent Tube is it still good?

- a) 30 minutes
- b) The mixture should be used immediately
- c) 24 hours
- d) 8 hours

7. When testing a nasopharyngeal swab specimen, what is the incubation time when the swab is in the Reagent solution?

- a) 1 minute
- b) 2 minutes
- c) There is no incubation time
- d) 10 minutes

8. What are the age limits of the QuickVue RSV 10 test?

- a) Less than 18 years of age
- b) There are none
- c) Less than 6 years of age
- d) Less than 6 months of age

9. What specimen types can be used when running this test?

- a) Nasopharyngeal swab only
- b) Nasopharyngeal swab, nasopharyngeal aspirate/wash
- c) Nasopharyngeal swab, nasal swab
- d) Nasopharyngeal swab, nasal swab, nasopharyngeal aspirate/wash

10. Which viral transport media is not compatible with this kit?

- a) M4-RT
- b) BD Universal Viral Transport Media
- c) Saline
- d) Copan Universal Transport Media

QuickVue® RSV 10 Quiz

Circle the correct answer

This quiz is an educational tool intended to assist facilities in evaluating their operators' understanding of the QuickVue RSV 10 test procedure. This quiz is not intended to be used as sole evidence of operator training or competency. Facilities are responsible for ensuring the quality of the testing performed by their operators. When testing controls or patient specimens, follow the current package insert instructions and/or procedure card included in the kit.

1. How often does the manufacturer recommend that external controls be run on the kit?

- a) Once per kit
- b) Once per each new shipment of kits
- c) Once per day
- d) Once per each untrained operator, once for each new shipment of kits – provided that each different lot in a shipment is tested

2. At what temperature do the kits need to be stored?

- a) -20°C (frozen)
- b) 15-30°C (room temp.)
- c) 2-30°C (refrigerated or room temp.)
- d) 2-8°C (refrigerated)

3. To guarantee the performance characteristics indicated in the package insert, which swab should be used when collecting a nasopharyngeal swab specimen?

- a) Cotton-tipped swabs
- b) Any type of swabs
- c) Any nasal swab
- d) The swabs provided in the kit

4. What is the read time for this kit?

- a) 10 minutes
- b) 15 minutes
- c) 20 minutes
- d) 8 minutes

5. Can the test be read after the read time?

- a) Yes, any time after the read time
- b) Yes, up to 5 minutes after the read time
- c) No
- d) Yes, up to 10 minutes after the read time

6. How long after adding the Reagent solution to the Reagent Tube is it still good?

- a) 30 minutes
- b) The mixture should be used immediately
- c) 24 hours
- d) 8 hours

7. When testing a nasopharyngeal swab specimen, what is the incubation time when the swab is in the Reagent solution?

- a) 1 minute
- b) 2 minutes
- c) There is no incubation time
- d) 10 minutes

8. What are the age limits of the QuickVue RSV 10 test?

- a) Less than 18 years of age
- b) There are none
- c) Less than 6 years of age
- d) Less than 6 months of age

9. What specimen types can be used when running this test?

- a) Nasopharyngeal swab only
- b) Nasopharyngeal swab, nasopharyngeal aspirate/wash
- c) Nasopharyngeal swab, nasal swab
- d) Nasopharyngeal swab, nasal swab, nasopharyngeal aspirate/wash

10. Which viral transport media is not compatible with this kit?

- a) M4-RT
- b) BD Universal Viral Transport Media
- c) Saline
- d) Copan Universal Transport Media