



1055 East State Street, Suite 100
Athens, OHIO 45701 USA
Phone: +1(740) 589-3300
Fax: +1(740) 592-9820
www.dhiusa.com

HSV Training Panel

HSV Panel, 24 samples (1-mL fill volume)

REF: 02-480024.v2

For *Training Use*

Please contact Diagnostic HYBRIDS Technical Services for technical assistance regarding this procedure.

USA Toll-free: (866) 344-3477
Direct: +1(740) 589-3300
Fax: +1(740) 592-9820
E-mail: technical_services@dhiusa.com

Symbols Lexicon/Glossary

- Temperature limit
- Consult instructions for use
- Manufacturer
- Authorized representative in the European Community
- Batch code/lot number
- Catalog number
- Use by YYYY-MON-DD
- CE mark of conformity (Conformité Européen)

Intended Use:

For Training use.

The *Herpes simplex Virus (HSV)* panel is comprised of 24 frozen 1-mL liquid aliquots. Six of these aliquots have been tested and confirmed to contain viable clinical isolates of *Herpes simplex virus* type 1. An additional six have been tested and confirmed to contain viable clinical isolates of *Herpes simplex virus* type 2. The remaining twelve aliquots have been tested and confirmed as negative for the presence of *Herpes simplex virus*. The panel has been designed for use with cell culture. The validation of this panel with non-cell culture systems is the responsibility of the end-user.

Principle:

Twelve HSV viruses have been propagated in cell culture, diluted and frozen at -70° C. The actual titer in each aliquot is unknown, and has not been quantitatively assayed. Inoculation of this panel should be performed as described in the Product Insert or Instructions for Use of the system in which it will be used.

Storage and handling instructions:

It is recommended that the panel be stored at -70°C or lower. Once the panel has been thawed it should not be re-frozen.

Warnings and Precautions:

- Personnel working with this panel must be properly trained in virus culture and safe handling techniques.¹
- Manipulations which present potential personnel hazards should be conducted in a Class II hood with gloves worn at all times and in accordance with the guidelines presented in the CDC-NIH Manual.

¹ Biosafety in Microbiological and Biomedical Laboratories (BMBL), 5th edition, 2007, CDC-NIH manual.
<http://www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm>

Virus Handling and Usage Instructions:

- The HSV panel should be transferred from the shipping container to a -70°C or lower freezer immediately upon receipt.
- At the time of use, rapidly thaw the panel aliquots in a 35° to 37°C water bath.
- Vortex the aliquots for 10-seconds prior to use.
- Inoculate the panel aliquots without further dilutions. (200-µL is the recommended volume for inoculation)

HSV Training Panel Key

Panel Member	Expected Result	Panel Member	Expected Result
1	Negative	13	Negative
2	HSV-2	14	HSV-2
3	HSV-1	15	Negative
4	Negative	16	Negative
5	Negative	17	HSV-2
6	HSV-2	18	Negative
7	Negative	19	HSV-1
8	HSV-1	20	Negative
9	HSV-2	21	HSV-1
10	HSV-1	22	Negative
11	Negative	23	Negative
12	HSV-2	24	HSV-1

Warranty Statement

These products are warranted to perform as described in their labeling and the Diagnostic HYBRIDS literature when used in accordance with their instructions. THERE ARE NO WARRANTIES WHICH EXTEND BEYOND THIS EXPRESS WARRANTY AND DIAGNOSTIC HYBRIDS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. Diagnostic HYBRIDS sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of Diagnostic HYBRIDS to repair or replace the products.

