

QuickVue® Influenza Nasal Swab Transport Tubes

25 foam-tipped nasal swabs inside plastic transport tubes.
For use with QuickVue Influenza and QuickVue Influenza A+B tests.
Cat. #20171



Specimens should be tested as soon as possible after collection. Samples may be stored refrigerated (2–8°C), or at room temperature (15–30°C), in a clean, dry, closed container (such as the tubes pictured above) for up to eight hours prior to testing.¹

Note: If dual testing for both influenza and RSV using the QuickVue Influenza A+B and RSV test kits, nasal swabs should not be used. Nylon flocked nasopharyngeal swabs made by Copan Diagnostics are recommended for use. Nasopharyngeal (NP) swabs may be ordered through various distributors, such as Fisher Scientific, Cardinal Health, Hardy Diagnostics, Diagnostic Hybrids, Inc. and Infolab, using Item #501CS01.US (NP swabs only) or Item #551C (NP swabs in transport tubes).

¹QuickVue Influenza package insert.



Worldwide Headquarters
10165 McKellar Court
San Diego, CA 92121 USA
www.quidel.com

To order: Call your Quidel Account Manager or your local distributor

QuickVue® Swab Transport Tube Set
Catalog #20171

For more information on QuickVue products, visit our website at www.quidel.com
or call your Quidel Account Manager:
800.874.1517 (U.S.) • 858.552-1100 (outside U.S.)