

Technical Bulletin

Synagis[®] does not interfere with the QuickVue[®] Influenza and RSV test kits

Synagis was evaluated with the QuickVue Influenza, QuickVue Influenza A+B, and QuickVue RSV test kits at a concentration of 1.0 mg/mL and did not interfere with test performance.

Please call Quidel Technical Support at 800-874-1517 (USA only), 858-552-1100 or technicalsupport@quidel.com if you have any questions regarding the QuickVue Influenza or RSV tests or any other Quidel product. Our hours of operation are Monday-Friday, 7:00 a.m.-5:00 p.m. Pacific Time.

You may also visit our website at www.quidel.com or www.flutest.com information on Quidel's line of Rapid Diagnostics, Bone Health and Autoimmune & Complement product lines. Other product information available on our website includes: CPT codes, CLSI, MSDS and package inserts.