

**Section 1 – Product and Company Identification****1.1 Manufacturer Information**

Quidel Corporation	Phone:	1.800.874.1517	Web:	<a href="http://quidel.com">quidel.com</a>
2005 East State Street, Suite 100	Fax:	1.740.592.9820	E-mail:	<a href="mailto:gehs@quidel.com">gehs@quidel.com</a>
Athens, OH 45701	Emergency (24-Hour):	1.866.519.4752		

**1.2 Product Information**

Product Name: Mouse PTH Antibody (53-84)  
(Catalog # 20-2310 for custom order quantities and Catalog # 21-2310 for fixed quantity of 100 µg)

Intended Use: Affinity Purified Goat Anti-Mouse PTH Antibody This component is for professional use only. Not intended for human or animal diagnostic or therapeutic uses.

**Section 2 – Hazard Identification**

**2.1 Classification of the Substance of Mixture** Not a hazardous substance or mixture.

**2.2 GHS Label Elements, including Precautionary Statements** Not a hazardous substance or mixture.

**2.3 Hazards Not Otherwise Classified (HNOC) or not covered by GHS**

Sodium azide (<0.1%) is present as a preservative. Avoid disposal of the product down sanitary or industrial plumbing systems. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides.

**Section 3 – Composition / Information on Ingredients****3.1 Mixtures**

- Contains non-hazardous quantities of proprietary ingredients according to OSHA (29 CFR1910.1200)
- Not a hazardous mixture according to Regulation (EC) No. 1272/2008.

**3.2** No chemicals need to be disclosed according to the applicable regulations for the components in this kit.

**Section 4 – First Aid Measures****4.1 Description of First Aid Measures**

**If inhaled:** Move the person to fresh air and support breathing as required.

**In case of skin contact:** Wash affected area with soap and water. Seek medical advice if irritation develops.

**In case of eye contact:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, get medical attention.

**If swallowed:** Never give anything by mouth to an unconscious person. Rinse mouth with water (do not swallow). Seek medical advice if irritation develops.

**4.2 Most Important Symptoms and Effects (both acute and delayed)**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**4.3 Indication of Any Immediate Medical Attention and Special Treatment Needed**

No specific measures identified.

**Section 5 – Fire Fighting Measures**

Only individuals properly trained and issued appropriate personal protective equipment should respond and attempt to extinguish a fire.

**Section 5 – Fire Fighting Measures (cont'd)**

- 5.1 Suitable Extinguishing Media**  
For small fires, use dry chemical, carbon dioxide or alcohol-resistant foam.
- 5.2 Special Hazards Arising from the Substance or Mixture**  
Nature of decomposition products not known.
- 5.3 General Fire Hazards**  
The components within this kit will not significantly contribute to the intensity of a fire.
- 5.4 Fire Fighting Equipment**  
Firefighters should wear full protective gear when responding to fires.

**Section 6 – Accidental Release Measures**

Only individuals properly trained and issued appropriate personal protective equipment should respond and attempt to clean up a spill or release. A large spill of the components contained within this kit is unlikely.

- 6.1 Personal Precautions** Avoid breathing vapors, mist or gas.
- 6.2 Environmental Precautions** Contain spill to prevent migration to drains, sewers or open water sources. Discharge to the environment must be avoided.
- 6.3 Methods and Materials for Clean-Up** Wipe up with inert absorbent material (e.g., paper towel, etc.). Thoroughly wash the area with soap and water after a spill or release clean-up.
- 6.4 Recovery and Neutralization** Collect spilled material and clean-up supplies and place in a sealed container for disposal. Refer to Section 13 for disposal guidance.

**Section 7 – Handling and Storage**

- 7.1 Specific Use** For professional use only – Not for use by the general public.
- 7.2 Precautions for Safe Handling**  
As with all chemical and biological substances, avoid getting the components within this kit ON YOU or IN YOU. Wash exposed areas thoroughly after using this kit. Do not eat or drink while using this kit. This kit should be handled only by qualified clinical or laboratory employees trained on the use of this kit and who are familiar with the potential hazards. Universal Precautions should be followed when handling and working with this kit. Keep out of reach of the general public.
- 7.3 Universal Precautions**  
No known test method can offer complete assurance that product(s) derived from human blood, serum or plasma used in components of this kit will not transmit infectious agents. All blood products should be treated as potentially infectious; use universal precautions. Users should ensure the use of proper personal protective equipment (PPE), covering of any cuts or abrasions on the skin (hands), and follow standard protocols for the decontamination of work surfaces when testing has been completed. Always wash hands thoroughly after handling the components within this kit.
- 7.4 Conditions for Safe Storage (including any incompatibilities)**  
To maintain efficacy, store according to the Package Insert instructions.
- 7.5 Incompatibilities**  
To maintain efficacy, store according to the Package Insert instructions.

**Section 8 – Exposure Controls and Personal Protection**

- 8.1 Exposure Limits** No data available for the components within this kit.
- 8.2 Exposure Controls**
- |                               |   |
|-------------------------------|---|
| Engineering Measures          | Use with adequate ventilation.              |
| Personal Protective Equipment |   |
| <i>Respiratory Protection</i> | None needed under normal conditions of use. |

**Section 8 – Exposure Controls and Personal Protection (cont'd)**

<i>Skin Protection</i>	Handle with appropriately rated chemical resistant gloves. Gloves should be inspected prior to use. Use proper glove technique to remove gloves to avoid contact with skin. Wash hands after handling the components within this kit.
<i>Eye Protection</i>	Wear safety glasses with side shields or goggles to prevent eye contact.
<i>Body Protection</i>	Use body protection appropriate for the task. A laboratory coat is recommended.
<i>Hygiene Measures</i>	Wash hands before use, after use and at the end of the workday.

**8.3 Environmental Exposure Controls** No special environmental controls are required.

**8.4 Special Notes** No data available.

**Section 9 – Physical and Chemical Properties**

Characteristic	Mouse PTH Antibody (53-84)
Boiling Point (°C)	No data available
Melting Point (°C)	No data available
Specific Gravity (H <sub>2</sub> O = 1)	No data available
Vapor Pressure (mm Hg)	No data available
Vapor Density (Air = 1)	No data available
Evaporation Rate (Ether = 1)	No data available
pH	7.0 to 7.5
Solubility in Water	No data available
Appearance and Odor	No data available

**Section 10 – Stability and Reactivity**

Characteristic	Mouse PTH Antibody (53-84)
Component Stability	No data available
Hazard Reaction Potential	No data available
Conditions to Avoid	No data available
Materials to Avoid	No data available
Hazardous Decomposition Products	No data available

**Section 11 – Toxicological Properties**
**11.1 Information on Toxicological Effects**

Acute Toxicity	No data available
Skin Corrosion / Irritation	No data available
Respiratory or Skin Sensitization	No data available
Generative Cell Mutagenicity	No data available
Carcinogenicity	No component of this kit present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by the ACGIH, IARC, NTP or OSHA.
Reproductive Toxicity	No data available
Teratogenicity	No data available

**Section 11 – Toxicological Properties (cont'd)**

## Specified Target Organ Toxicity

Single Exposure: No data available

Repeated Exposure: No data available

Aspiration Hazard No data available

## Potential Health Effects

*Inhalation:* May cause respiratory tract irritation.

*Skin:* May cause skin irritation upon contact.

*Ingestion:* May be harmful if swallowed.

*Eyes:* May cause eye irritation.

**11.2 Signs and Symptoms of Exposure**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated for the components within this kit.

**11.3 Additional Information**

Liver – Irregularities – Based on human evidence (sodium azide).

**Section 12 – Ecological Information**
**12.1 Toxicity**

No data available

**12.2 Persistence and Degradability**

No data available

**12.3 Bioaccumulative Potential**

No data available

**12.4 Mobility in Soil**

No data available

**12.5 PBT and vPvB Assessment**

No data available

**12.6 Other Adverse Effects**

An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

**Section 13 – Disposal Considerations**
**13.1 Waste Disposal Instructions**

Utilize appropriate personal protective equipment and spill control when handling wastes generated from using this kit. Do not discharge any of the solutions, reagents or controls into drains, water courses or onto the ground.

**13.2 Disposal of Product and Contaminated Packaging**

Dispose of waste materials, unused components and contaminated packaging in compliance with country, federal, state and local regulations. If unsure of the applicable regulatory requirements, contact a licensed professional waste disposal service to dispose of this material.

**Section 14 – Transportation Information**
**14.1 U.S Department of Transportation (DOT)**

This component is not regulated for transport.

**14.2 International Air Transportation Association (IATA)**

This component is not regulated for transport.

**14.3 International Maritime Dangerous Goods (IMDG)**

This component is not regulated for transport.

**Section 15 – Regulatory Information**
**15.1 U.S. Federal Regulations**

OSHA Hazards None

SARA 302 The following chemicals are subject to reporting levels established by Sara Title III, Section 302:

Sodium Azide CAS #: 26628-22-8 Revision Date: 2007-07-01

SARA 313 The following chemicals are subject to reporting levels established by SARA Title III, Section 313:

Sodium Azide CAS #: 26628-22-8 Revision Date: 2007-07-01

SARA311/312 Hazards Sodium Azide CAS #: 26628-22-8 Acute Health Hazard; Chronic Health Hazard

**Section 15 – Regulatory Information (cont'd)**
**15.2 State Regulations**

The following chemicals appear on one or more of the following state hazardous substance lists:

Chemical Name	CAS #	CA	MA	MN	NJ	PA	RI
Sodium Azide	26628-22-8	Yes	Yes	Yes	Yes	Yes	Yes

California Prop 65: This kit does not contain any chemicals known to the State of California to cause cancer, birth defects, or any other reproductive harm.

**15.3 Canadian – WHMIS Classification**

Chemical Name	CAS #	Classification
--	--	--

**15.4 Additional Regulatory Information**
**Safety, Health and Environmental Regulations/Legislation Specific for the Mixture**


European Union – Regulation (EC) No. 1907/2006

**Chemical Safety Assessment**

Not completed for the components contained within this kit.

**NFPA Kit Classification:**

Health Hazard (blue): 1  
 Fire Hazard (red): 0  
 Reactivity (yellow): 0  
 Special Hazards (white): None


**Section 16 – Other Information**

**PREPARED BY:** Quidel Corporation  
 12544 High Bluff Drive, Suite 200  
 San Diego, CA 92130

**SUPERSEDES:** November 21, 2017

**REVISIONS:** Updated product name and document number.

The information above is provided in good faith. It is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability, fitness for a particular purpose or of any other type, expressed or implied, with respect to products described or data or information provided, and we assume no liability resulting from the use of such products, data or information. Users should make their own investigations to determine the suitability of the information for their particular purposes, and the user assumes all risk arising from their use of the material. The user is required to comply with all laws and regulations relating to the purchase, use, storage and disposal of the material, and must be familiar with and follow generally accepted safe handling procedures. In no event shall Quidel be liable for any claims, losses, or damages of any individual or for lost profits or any special, indirect, incidental, consequential or exemplary damages of any kind, howsoever arising, even if Quidel has been advised of the possibility of such damages.

SD2310000EN00 (08/18)