



A replacement component for Quidel Product 01-010000.v2 D³ Ultra™ DFA Respiratory Virus Screening & ID Kit.

FOR *IN VITRO* DIAGNOSTIC USE



INTENDED USE

The Mounting Fluid is a replacement component for Quidel Product 01-010000.v2 D³ Ultra DFA Respiratory Virus Screening & ID Kit.

MATERIALS PROVIDED

Quidel Product 01-002015 contains:

Mounting Fluid

15 mL

One dropper bottle containing 15 mL of an aqueous, buffered, stabilized solution of glycerol and 0.1% sodium azide.

STORAGE AND HANDLING

Reagent will retain its full potency through the expiration date when stored at 2°C to 8°C. Store in the dark.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use
- Use undiluted.
- Refer to 01-010000.v2 Package Insert for complete directions.
- Sodium azide is included in this reagent at 0.1%.
 - ▶ Aqueous solutions of sodium azide, when mixed with acids, may liberate toxic gas.
 - ▶ Sodium azide may react with lead and copper plumbing to form highly explosive metal azides.
 - ▶ Avoid disposal of these solutions down sanitary or industrial plumbing systems.
 - ▶ Avoid release to the environment.
- Testing should be performed in an area with adequate ventilation.
- Dispose of containers and unused contents in accordance with Federal, State and Local regulatory requirements.
- Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this kit.
- Wash hands thoroughly after handling.
- For additional information on hazard symbols, safety, handling and disposal of the components within this kit, please refer to the Safety Data Sheet (SDS) located at quidel.com.

ASSISTANCE

To place an order or for technical support, please contact a Quidel Representative at 800.874.1517 (in the U.S.) or 858.552.1100 (outside the U.S.), Monday through Friday, from 8:00 a.m. to 5:00 p.m., Eastern Time. Orders may also be placed by fax at 740.592.9820. For e-mail support contact customerservice@quidel.com or technicalsupport@quidel.com.

For services outside the U.S.A., please contact your local distributor. Additional information about Quidel, our products, and our distributors can be found on our website quidel.com.

REF

01-002015 – Mounting Fluid

IVD



EC REP

MDSS GmbH
Schiffgraben 41
30175 Hannover,
Germany



Diagnostic Hybrids, Inc. – a subsidiary of Quidel Corporation
2005 East State Street, Suite 100
Athens, OH 45701 USA
quidel.com

PI167MF00EN00 (10/18)

GLOSSARY

REF

Catalogue number



CE mark of conformity

EC REP

Authorized Representative
in the European Community

LOT

Batch code



Use by



Manufacturer



Temperature limitation



Intended use



Consult e-labeling
instructions for use



Do not reuse

IVD

For *In Vitro* diagnostic use
