



QUIDEL

Pen/Strep 100X



INTENDED USE

One bottle containing 100 mL of a 0.9% sodium chloride solution with 10,000 units/mL of Penicillin and 10,000 units/mL of Streptomycin for use in the preparation of cell culture media. This solution has been tested for the presence of endotoxins as well as for its compatibility with cell cultures.

MATERIALS PROVIDED

Cat. #10-010101

Component	Quantity	Storage
Pen/Strep 100X	1 bottle (100 mL)	2°C to 8°C

STORAGE/HANDLING

Reagent will retain its full potency through the expiration date when stored at 2°C to 8°C. Store in the dark.

INSTRUCTIONS

Dilute to 1X before use.

WARNINGS AND PRECAUTIONS

- General purpose laboratory reagent. Not intended for human or animal diagnostic or therapeutic uses
- May be an eye, skin or respiratory irritant. For eye exposure, rinse with water for several minutes. For skin exposure, wash with soap and water. For respiratory difficulties, remove to fresh air. Seek medical attention if any condition persists.
- Testing should be performed in an area with adequate ventilation.
- Dispose of containers and unused contents in accordance with Federal, State and Local regulatory requirements.
- Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this kit.
- Wash hands thoroughly after handling.
- For additional information on hazard symbols, safety, handling and disposal of the components within this kit, please refer to the Safety Data Sheet (SDS) located at quidel.com.

CUSTOMER AND TECHNICAL SUPPORT

To place an order or for Technical Support, please contact a Quidel Representative at 800.874.1517 (in the U.S.) or 858.552.1100 (outside the U.S.), Monday through Friday, from 8:00 a.m. to 5:00 p.m., Eastern Time. Orders may also be placed by fax at 740.592.9820. For e-mail support contact customerservice@quidel.com or technicalsupport@quidel.com.

For services outside the U.S., please contact your local distributor. Additional information about Quidel, our products, and our distributors can be found on our website quidel.com.

REF

10-010101 – Pen/Strep 100X

IVD



Diagnostic Hybrids, Inc. – a subsidiary of Quidel Corporation

2005 East State Street, Suite 100

Athens, OH 45701 USA

quidel.com

PI9010000EN00 (12/18)

GLOSSARY

REF

Catalogue number

LOT

Batch code



Use by



Manufacturer



Temperature limitation



Intended use



Consult e-labeling instructions for use



Do not reuse

IVD

For *in vitro* diagnostic use

CONT

Contents/Contains
