

### SECTION 1: Identification

#### 1.1. Identification

Product form : Mixture  
 Product name : Monoclonal Antibodies - A800 Series  
 Product code(s)

Product Description	Part Number	Quantity in Each Unit	Bulk Order Part Number	Quantity in Each Unit
Rat mAb Anti-Mouse Factor B (Intact) 5C1	A800	0.1 to 1.5 ml	IR049.192	2.0 to 50 ml
Rat mAb Anti-Mouse Factor B (Intact/Ba) 7D10	A801	0.1 to 1.5 ml	IR049.194	2.0 to 50 ml
Rat mAb Anti-Mouse Factor B (Intact/Bb) 5E10	A802	0.1 to 1.5 ml	IR049.193	2.0 to 50 ml
Rat mAb Anti-Mouse C3 (Intact/C3a) 8D11	A803	0.1 to 1.5 ml	IR049.180	2.0 to 50 ml
Rat mAb Anti-Mouse C3a (Neoepitope) 8G9	A804	0.1 to 1.5 ml	IR049.181	2.0 to 50 ml

#### 1.2. Recommended use and restrictions on use

Recommended use : For research use only, For medical diagnostic use.  
 Restrictions on use : Restricted to professional users, Professional Use of Medical Devices

#### 1.3. Supplier

##### Manufacturer

Quidel Corporation  
 2005 East State Street, Suite 100  
 Athens, 45701 - USA  
 T 1.800.874.1517 - F 1.740.592.9820  
[gehs@quidelortho.com](mailto:gehs@quidelortho.com) - [quidel.com](http://quidel.com)

#### 1.4. Emergency telephone number

Emergency number : 1.866.519.4752

### SECTION 2: Hazard(s) identification

#### 2.1. Classification of the substance or mixture

##### GHS US classification

Not classified

#### 2.2. GHS Label elements, including precautionary statements

According the corresponding national regulations there is no labelling obligation for this product.

#### 2.3. Other hazards which do not result in classification

No additional information available

#### 2.4. Unknown acute toxicity (GHS US)

Not applicable

### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	GHS US classification
sodium azide	(CAS No) 26628-22-8	0.02	Acute Tox. 2 (Oral), H300 Acute Tox. 1 (Dermal), H310 Aquatic Acute 1, H400 Aquatic Chronic 1, H410

Full text of hazard classes and H-statements : see section 16

### SECTION 4: First-aid measures

#### 4.1. Description of first aid measures

- First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing.
- First-aid measures after skin contact : Wash skin with plenty of water.
- First-aid measures after eye contact : Rinse eyes with water as a precaution.
- First-aid measures after ingestion : Call a poison center or a doctor if you feel unwell.

#### 4.2. Most important symptoms and effects (acute and delayed)

No additional information available

#### 4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

### SECTION 5: Fire-fighting measures

#### 5.1. Suitable (and unsuitable) extinguishing media

- Suitable extinguishing media : Water spray. Dry powder. Foam. Carbon dioxide.

#### 5.2. Specific hazards arising from the chemical

- Reactivity : The product is non-reactive under normal conditions of use, storage and transport.

#### 5.3. Special protective equipment and precautions for fire-fighters

- Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### 6.1.1. For non-emergency personnel

- Emergency procedures : Ventilate spillage area.

##### 6.1.2. For emergency responders

- Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

#### 6.2. Environmental precautions

Avoid release to the environment.

#### 6.3. Methods and material for containment and cleaning up

- Methods for cleaning up : Take up liquid spill into absorbent material.
- Other information : Dispose of materials or solid residues at an authorized site.

#### 6.4. Reference to other sections

For further information refer to section 13.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

- Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment.
- Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

#### 7.2. Conditions for safe storage, including any incompatibilities

- Storage conditions : Store in a well-ventilated place. Keep cool.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

sodium azide (26628-22-8)		
ACGIH	ACGIH OEL C	0.29 mg/m <sup>3</sup>
ACGIH	ACGIH OEL C [ppm]	0.11 ppm

### 8.2. Appropriate engineering controls

- Appropriate engineering controls : Ensure good ventilation of the work station.  
 Environmental exposure controls : Avoid release to the environment.

### 8.3. Individual protection measures/Personal protective equipment

#### Materials for protective clothing:

Lab coat

#### Hand protection:

Wear protective gloves. Wash your hands

#### Eye protection:

Safety glasses

#### Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment

#### Personal protective equipment symbol(s):



#### Other information:

Do not eat, drink or smoke when using this product.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Clear, colorless liquid
Color	: Clear, colorless
Odor	: Odorless
Odour threshold	: No data available
pH	: 7.2
Melting point	: Not applicable
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Flammability	: Not applicable.
Vapour pressure	: No data available
Relative vapour density at 20°C	: No data available
Relative density	: No data available
Solubility	: No data available
Partition coefficient n-octanol/water (Log Pow)	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive limits	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

### 10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified  
 Acute toxicity (dermal) : Not classified  
 Acute toxicity (inhalation) : Not classified

<b>sodium azide (26628-22-8)</b>	
LD50 oral rat	27 mg/kg bodyweight (Rat, Experimental value, Oral)
LD50 dermal rabbit	19 – 48 mg/kg bodyweight (Rabbit, Inconclusive, insufficient data, Dermal)
LC50 Inhalation - Rat	0.054 – 0.52 mg/l air Animal: rat, Guideline: EPA OPPTS 870.1300 (Acute inhalation toxicity)
ATE US (oral)	27 mg/kg bodyweight
ATE US (dermal)	19 mg/kg bodyweight

Skin corrosion/irritation : Not classified  
 pH: 7.2

Serious eye damage/irritation : Not classified  
 pH: 7.2

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : Not classified

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Viscosity, kinematic : No data available

### SECTION 12: Ecological information

#### 12.1. Toxicity

Ecology - general : The product is not considered harmful to aquatic organisms nor to cause long-term adverse effects in the environment.

sodium azide (26628-22-8)	
LC50 - Fish [1]	2.75 – 3.28 mg/l (Equivalent or similar to OECD 203, 96 h, Oncorhynchus mykiss, Flow-through system, Fresh water, Experimental value)
LC50 - Other aquatic organisms [1]	1 (1 – 10) mg/l (96 h)
EC50 - Crustacea [1]	4.2 mg/l Source: NCIS
EC50 - Other aquatic organisms [1]	5 (5 – 14) mg/l (Protozoa; TOXICITY TEST)

#### 12.2. Persistence and degradability

sodium azide (26628-22-8)	
Persistence and degradability	Biodegradability: not applicable.
Biochemical oxygen demand (BOD)	Not applicable
Chemical oxygen demand (COD)	Not applicable
ThOD	Not applicable
BOD (% of ThOD)	Not applicable

#### 12.3. Bioaccumulative potential

sodium azide (26628-22-8)	
Partition coefficient n-octanol/water (Log Pow)	0.16 Source: NIOSH
Bioaccumulative potential	not bioaccumulable.

#### 12.4. Mobility in soil

sodium azide (26628-22-8)	
Surface tension	No data available (test not performed)
Organic Carbon Normalized Adsorption Coefficient (Log Koc)	2.63 (log Koc, Calculated value)
Ecology - soil	Low potential for adsorption in soil.

#### 12.5. Other adverse effects

No additional information available

### SECTION 13: Disposal considerations

#### 13.1. Disposal methods

Regional legislation (waste) : Disposal must be done according to official regulations.  
 Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.  
 Sewage disposal recommendations : Disposal must be done according to official regulations.  
 Ecology - waste materials : Avoid release to the environment.

### SECTION 14: Transport information

#### Department of Transportation (DOT)

In accordance with DOT

Not regulated

#### Transport by sea

Not regulated

#### Air transport

Not regulated

### SECTION 15: Regulatory information

#### 15.1. US Federal regulations

<b>sodium azide (26628-22-8)</b>	
Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	1000 lb
RQ (Reportable quantity, section 304 of EPA's List of Lists)	1000 lb
SARA Section 302 Threshold Planning Quantity (TPQ)	500 lb

#### 15.2. International regulations

##### CANADA

<b>sodium azide (26628-22-8)</b>
Listed on the Canadian DSL (Domestic Substances List)

##### EU-Regulations

No additional information available

##### National regulations

<b>sodium azide (26628-22-8)</b>
Listed on the United States TSCA (Toxic Substances Control Act) inventory - Status: Active
Listed on INSQ (Mexican National Inventory of Chemical Substances)

#### 15.3. US State regulations

California Proposition 65 - This product does not contain any substances known to the state of California to cause cancer, developmental and/or reproductive harm

### SECTION 16: Other information

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Revision date : 04/18/2023

Other information : Please review product insert prior to using this product.

Full text of H-statements:

H300	Fatal if swallowed.
H310	Fatal in contact with skin.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.

Indication of changes:

Added.

SDS US Quidel

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.*