

### Section 1 – Product and Company Identification

**1.1 Manufacturer Information**

Quidel Corporation	Phone:	1.800.874.1517	Web:	<a href="http://guidel.com">guidel.com</a>
2005 East State Street, Suite 100	Fax:	1.740.592.9820	E-mail:	<a href="mailto:gehs@guidel.com">gehs@guidel.com</a>
Athens, OH 45701	Emergency (24-Hour):	1.866.519.4752		

**1.2 Product Information**

Product Name: Solana Influenza A+B test kit

 Intended Use: An *in vitro* diagnostic test, using isothermal amplification technology (helicase-dependent amplification, HDA), for the qualitative detection and differentiation of influenza A and influenza B viral RNA in nasal and nasopharyngeal swabs from patients with signs and symptoms of respiratory infection. Intended for use only with the Solana® instrument.

Components: Process Buffer (1.55 mL), Reaction tubes (Lyophilized cake in 0.2 mL tube)

### Section 2 – Hazard Identification

- |  |                                       |
|--|---------------------------------------|
| <b>2.1 Classification of the Substance of Mixture(s)</b>                 | Not a hazardous substance or mixture. |
| <b>2.2 GHS Label Elements, including Precautionary Statements</b>        | Not a hazardous substance or mixture. |
| <b>2.3 Hazards Not Otherwise Classified (HNOC) or not covered by GHS</b> | None                                  |

### Section 3 – Composition / Information On Ingredients

**3.1 Mixtures**

Contains non-hazardous quantities of proprietary ingredients according to OSHA (29 CFR 1910.1200).

Not a hazardous mixture according to Regulation (EC) No. 1272/2008.

Component	Chemical Name	CAS #	EINECS #	Concentration (%)
Process Buffer	Glycerol	56-81-5	200-289-5	<5.5

### Section 4 – First Aid Measures

**4.1 Description of First Aid Measures**

**If inhaled:** Move the person to fresh air and support breathing as required.

**In case of skin contact:** Wash affected area with soap and water. Seek medical advice if irritation develops.

**In case of eye contact:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists get medical attention.

**If swallowed:** Never give anything by mouth to an unconscious person. Rinse mouth with water. Seek medical advice if irritation develops.

**4.2 Most Important Symptoms and Effects (both acute and delayed)**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**4.3 Indication of Any Immediate Medical Attention and Special Treatment Needed**

No specific measures identified.

### Section 5 – Fire Fighting Measures

Only individuals properly trained and issued appropriate personal protective equipment should respond and attempt to extinguish a fire.

**5.1 Suitable Extinguishing Media**

For small fires, use dry chemical, carbon dioxide or alcohol-resistant foam.

**5.2 Special Hazards Arising From the Substance or Mixture**

Nature of decomposition products not known.

**5.3 General Fire Hazards**

The components within this kit will not significantly contribute to the intensity of a fire.

**5.4 Fire Fighting Equipment**

Firefighters should wear full protective gear when responding to fires.

### Section 6 – Accidental Release Measures

Only individuals properly trained and issued appropriate personal protective equipment should respond and attempt to clean up a spill or release. A large spill of the components contained within this kit is unlikely.

**6.1 Personal Precautions**

Avoid breathing vapors, mist or gas.

**6.2 Environmental Precautions**

Contain spill to prevent migration to drains, sewers or open water sources. Discharge to the environment must be avoided.

**6.3 Methods and Materials for Clean-Up**

Soak up with inert absorbent material (e.g., paper towel, etc.). Thoroughly wash the area with soap and water after a spill or release clean-up.

**6.4 Recovery and Neutralization**

Collect spilled material and clean-up supplies and place in a sealed container for disposal. Refer to Section 13 for disposal guidance.

### Section 7 – Handling and Storage

**7.1 Specific Use**

For professional use only – Not for use by the general public.

**7.2 Precautions for Safe Handling**

As with all chemical and biological substances, avoid getting the components within this kit ON YOU or IN YOU. Wash exposed areas thoroughly after using this kit. Do not eat or drink while working with this kit. This kit should be handled only by qualified clinical or laboratory employees trained on the use of these controls, reagents, solutions or buffers and who are familiar with the potential hazards. **Universal Precautions should be followed when handling and working with this kit.** Keep out of reach of the general public.

**7.3 Conditions for Safe Storage (including any incompatibilities)**

To maintain efficacy, store according to the package insert instructions.

**7.4 Incompatibilities**

To maintain efficacy, store according to the package insert instructions.

### Section 8 – Exposure Controls and Personal Protection

**8.1 Exposure Limits**

No data available for the components within this kit.

**8.2 Exposure Controls**

Engineering Measures

Use with adequate ventilation.

Personal Protective Equipment

*Respiratory Protection*

None needed under normal conditions of use.

*Skin Protection*

Handle with appropriately rated chemical resistant gloves. Gloves should be inspected prior to use. Use proper glove technique to remove gloves to avoid contact with skin. Wash hands after handling the components within this kit.

**Section 8 – Exposure Controls and Personal Protection (cont'd.)**

*Eye Protection*                      Wear safety glasses with side shields or goggles to prevent eye contact.  
*Body Protection*                    Use body protection appropriate for the task. A laboratory coat is recommended.  
*Hygiene Measures*                 Wash hands before use, after use and at the end of the workday.

**8.3 Environmental Exposure Controls**    No special environmental controls are required.

**8.4 Special Notes**                         No data available.

**Section 9 – Physical and Chemical Properties**

Characteristic	Process Buffer	Reaction Tubes
Boiling Point (°C)	No data available	No data available
Melting Point (°C)	No data available	No data available
Specific Gravity (H <sub>2</sub> O = 1)	No data available	No data available
Vapor Pressure (mm Hg)	No data available	No data available
Vapor Density (Air = 1)	No data available	No data available
Evaporation Rate (Ether = 1)	No data available	No data available
pH	8.0	No data available
Solubility in Water	Soluble	Soluble
Appearance and Odor	Clear, Odorless	White cake, Odorless

**Section 10 – Stability and Reactivity**

Characteristic	Process Buffer	Reaction Tubes
Component Stability	Stable	Stable
Hazard Reaction Potential	No data available	No data available
Conditions to Avoid	No data available	No data available
Materials to Avoid	No data available	No data available
Hazardous Decomposition Products	No data available	No data available

**Section 11 – Toxicological Properties**
**11.1 Information on Toxicological Effects**

Acute Toxicity                                No data available  
Skin Corrosion / Irritation                No data available  
Respiratory or Skin Sensitization        No data available  
Generative Cell Mutagenicity             No data available  
Carcinogenicity                                No component of this kit present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by the ACGIH, IARC, NTP or OSHA.  
Reproductive Toxicity                        No data available  
Teratogenicity                                 No data available

### Section 11 – Toxicological Properties (cont'd.)

**Specified Target Organ Toxicity**

Single Exposure: No data available      Repeated Exposure: No data available

Aspiration Hazard      No data available

**Potential Health Effects**

*Inhalation:*      May cause respiratory tract irritation.      *Skin:*      May cause skin irritation upon contact.  
*Ingestion:*      May be harmful if swallowed.      *Eyes:*      May cause eye irritation.

**11.2 Signs and Symptoms of Exposure**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated for the components within this kit.

**11.3 Additional Information**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated for the components within this kit.

Kidney – Irregularities – Based on Human Evidence (Glycerol)

### Section 12 – Ecological Information

<b>12.1 Toxicity</b>	No data available
<b>12.2 Persistence and Degradability</b>	No data available
<b>12.3 Bioaccumulative Potential</b>	No data available
<b>12.4 Mobility in Soil</b>	No data available
<b>12.5 PBT and vPvB Assessment</b>	Assessment not available as chemical safety assessment not conducted.
<b>12.6 Other Adverse Effects</b>	An environmental hazard cannot be excluded in the event of unprofessional handling or disposal. Avoid disposal to the environment.

### Section 13 – Disposal Considerations

**13.1 Waste Disposal Instructions**

Utilize appropriate personal protective equipment and spill control when handling wastes generated from using this kit. Do not discharge any of the controls, buffers, reagents or solutions into drains, water courses or onto the ground.

**13.2 Disposal of Product and Contaminated Packaging**

Dispose of waste materials, unused components and contaminated packaging in compliance with country, federal, state and local regulations. If unsure of the applicable regulatory requirements, contact a licensed professional waste disposal service to dispose of this material.

### Section 14 – Transportation Information

<b>14.1 U.S Department of Transportation (DOT)</b>	This kit not regulated for transport.
<b>14.2 International Air Transportation Association (IATA)</b>	This kit not regulated for transport.
<b>14.3 Canada</b>	This kit not regulated for transport.

### Section 15 – Regulatory Information

**15.1 U.S. Federal Regulations**

OSHA Hazards      None  
 SARA 302      The following chemicals are subject to reporting levels established by Sara Title III, Section 302:  
                          No chemicals in this kit are subject to the reporting requirements of SARA Title III, Section 302.

**Section 15 – Regulatory Information (cont'd.)**

SARA 313 The following chemicals are subject to reporting levels established by SARA Title III, Section 313:  
No chemicals in this kit are subject to the reporting requirements of SARA Title III, Section 313.

SARA311/312 Hazards No SARA hazards.

**15.2 State Regulations**

The following chemicals appear on one or more of the following state hazardous substance lists:

Chemical Name	CAS #	CA	MA	MN	NJ	PA	RI
Glycerol	56-81-5	No	Yes	Yes	Yes	Yes	Yes

California Prop 65: This product does not contain any chemicals known to the State of California to cause cancer, birth defects, or any other reproductive harm.

**15.3 Canadian – WHMIS**

Chemical Name	CAS #	Minimum Concentration
--	--	--

**15.4 Additional Regulatory Information**

**Safety, Health and Environmental Regulations/Legislation Specific for the Mixture** No data available

**Chemical Safety Assessment** Not completed for the components contained within this kit.

**HMIS Kit Classification:**

Health Hazard: 1  
Chronic Health Hazard: \*  
Flammability: 0  
Physical Hazard: 0

**NFPA Kit Classification:**

Health Hazard (blue): 1  
Fire Hazard (red): 0  
Reactivity (yellow): 0  
Special Hazards (white): None

*\*Use additional care when handling the components within this kit.*

**Section 16 – Other Information**

Every effort has been made to adhere to the hazard criteria and content requirements of the U.S. OSHA Hazard Communication Standard, European Communities Safety Data Sheets Directive, Canadian Controlled Products Regulations, UK Chemical Hazard information and Packaging Regulations, and UN Globally Harmonized System of Classification and Labeling of Chemicals.

**PREPARED BY:** Quidel Corporation  
12544 High Bluff Drive, Suite 200  
San Diego, CA 92130

**SUPERSEDES:** New

**REVISIONS:** Not applicable

The information above is provided in good faith. It is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability, fitness for a particular purpose or of any other type, expressed or implied, with respect to products described or data or information provided, and we assume no liability resulting from the use of such products, data or information. Users should make their own investigations to determine the suitability of the information for their particular purposes, and the user assumes all risk arising from their use of the material. The user is required to comply with all laws and regulations relating to the purchase, use, storage and disposal of the material, and must be familiar with and follow generally accepted safe handling procedures. In no event shall Quidel be liable for any claims, losses, or damages of any individual or for lost profits or any special, indirect, incidental, consequential or exemplary damages of any kind, howsoever arising, even if Quidel has been advised of the possibility of such damages.