



# QuickVue®

Results. Right. Now.

## Chlamydia



## Test and Treat Today®!

- Results in 12 minutes, or less.
- Most patients tested can be diagnosed and treated in just one visit .
- Simple extraction procedure with just two reagents.
- After extraction, just add the sample to the Test Cassette.
- All components necessary to perform test are contained in the kit, including swabs.
- Room temperature storage (9 months from date of manufacture).
- Built-in Positive and Negative controls confirm the test has performed properly.
- Positive and Negative controls provided in the kit.

# Procedure:

