

DFA Herpes Simplex Virus



Identification and typing kits

A Complete System Approach

- Combine the D³ DFA HSV Identification or Identification and Typing Kit with H&V-Mix MixedCells, FLOQSwabs and UTM for a complete system approach to your Herpes simplex virus testing needs.




Rapid Turnaround Time

- Most lesion or genital herpes cultures show positive results within 24 hours on H&V-Mix MixedCells. Confirm identification with one 15-minute staining step.

Flexible Reagent Formats

- Order as a complete kit containing reagents, control slides, PBS and mounting fluid, or order reagents independently as an Identification and Typing Duo Pack (HSV-1 and HSV-2 Reagents packaged together) or Single HSV Identification Reagent. It's your choice!

Recommended protocol – Shell Vial



H&V-Mix Products

Vials

98-0102 Shell Vial – With Coverslip

Plates

98-2412 24 well/12 fill
 98-2418 24 well/18 fill
 98-2424 24 well/24 fill
 98-4812 48 well/12 fill
 98-4818 48 well/18 fill
 98-4824 48 well/24 fill

Standard 2% FBS Refeed Medium

10-320100 100 mL
 10-320500 500 mL

Hs27 ReadyCells® System

Hs27ReadyCells

F-87-0102-24 24 Shell Vials – With Coverslip

ReadyCells Heatblock

99-610-009 ReadyCells Heatblock (120 volts)
 (used for rapid, controlled thawing of ReadyCells vials)

Individual Components

01-090001	D ³ DFA Herpes Simplex Virus Identification & Typing Reagent Duo Pack Includes: HSV-1 DFA Reagent (5 mL), HSV-2 DFA Reagent (5 mL)	20-4810020	HSV-1 Virus Control Stock (Single Use, 400 µL) (shipped frozen)
01-085005	D ³ DFA Herpes Simplex Virus Reagent (5 mL)	20-4820020	HSV-2 Virus Control Stock (Single Use, 400 µL) (shipped frozen)
01-00100	HSV-1/HSV-2 Antigen Control Slides (10 slides per set)	02-480024.v2	Virus Training HSV Panel (24 vials, 1 mL fill volume) (shipped frozen)
01-090025	40X PBS Concentrate (25 mL)	02-461024	H&V Training Panel (24 vials, 1 mL fill volume) (shipped frozen)
01-002007b	Mounting Fluid (7 mL)		

D³ DFA Herpes Simplex Virus Identification Kit – Cat. #01-080000

Includes: HSV DFA Reagent (5 mL), 40X PBS Concentrate (25 mL), Mounting Fluid (7 mL), HSV-1/HSV-2 Antigen Control Slides (5 slides)

D³ DFA Herpes Simplex Virus Identification & Typing Kit – Cat. #01-090000

Includes: HSV-1 DFA Reagent (5 mL), HSV-2 DFA Reagent (5 mL), 40X PBS Concentrate (25 mL), Mounting Fluid (7 mL), HSV-1/HSV-2 Antigen Control Slides (5 slides)

SL1102 | SS1013000EN00 (04/17)

