

IFA Enterovirus



Rapid Enterovirus testing solutions

Sensitive And Enterovirus-Specific Identification Stain

- Superior staining with proven proprietary monoclonal antibodies
- A blend of monoclonal antibodies to assure broad enterovirus coverage
- Antibodies targeting the VP1 protein allow for enterovirus-specific identification
- Clear staining results without a blocking reagent


Simplified Workflow

- No need to read CPE over several days
- Complete kit format simplifies Quality Control
- Super E-Mix combines BGMK-hDAF and A549 cells to allow isolation of enterovirus in a single culture system
- Super E-Mix now available in a frozen ReadyCell format – ready for use in 4 minutes with long storage capabilities


Rapid Final Results For All Patients

- Results in as little as 48 hours
- Definitive final results for enterovirus without additional confirmatory testing

Super E-Mix Cells infected with
Enterovirus 71, stained with
D³ IFA Enterovirus ID Kit



Recommended protocol



Individual Components

01-050001	D ³ Enterovirus Duo Pack Includes: Enterovirus MAb Reagent (5 mL), Anti-Mouse Conjugate (5 mL)
01-00080	Enterovirus Antigen Control Slides (10 Slides)
01-002007b	Mounting Fluid (7 mL) pH 8.0-8.4
01-057005	D ³ Enterovirus 71 DFA Reagent
01-058005	D ³ Enterovirus 71/Coxsackie A16 DFA Reagent

Super E-Mix

Vials	
92-0102	Shell Vial – With Coverslip
Plates	
92-2412	Cluster Plate, 24 well/12 fill
92-2418	Cluster Plate, 24 well/18 fill
92-2424	Cluster Plate, 24 well/24 fill
92-4812	Cluster Plate, 48 well/12 fill
92-4818	Cluster Plate, 48 well/18 fill
92-4824	Cluster Plate, 48 well/24 fill

Super E-Mix RefeedMedium
10-380100 100 mL

Enterovirus Detection Supplies

Virus Control Stocks

(Single Use, 400 µL) (shipped frozen)
20-4420010 Enterovirus-Coxsackie B
20-4510130 Enterovirus-Echovirus
02-440020 Enterovirus Training Panel, 24 samples
(1 mL fill volume) (shipped frozen)

Super E-Mix – ReadyCells®

Super E-Mix ReadyCells

F-92-0102-24 24 Shell Vials – With Coverslip

Standard 2% FBS Refeed Medium

10-320100 100 mL
10-320500 500 mL

ReadyCells Heatblock

99-610-009 ReadyCells Heatblock (120 volts)
(used for rapid, controlled thawing
of ReadyCells vials)

D³ IFA Enterovirus Identification Kit – Cat. #01-050000

Includes: Enterovirus MAb Reagent (5 mL), Anti-Mouse Conjugate (5 mL), 40X PBS Concentrate (25 mL),
Mounting Fluid (7 mL), Enterovirus Antigen Control Slides (5)

SL0701 | SS1780000EN00 (04/17)

