



QuickVue® TLI

Lactoferrin




Test and Treat Today®!


Aid in Diagnosis and Treatment of IBD

- Detects intestinal inflammation and helps to differentiate Irritable Bowel Syndrome (IBS) from active Inflammatory Bowel Disease (IBD).
- Determine elevated lactoferrin levels in patients and help guide treatment.¹
- Clear results in less than 10 minutes.
- Can be used with samples stored at room temperature for up to 2 weeks.
- Each kit comes with everything needed to perform the test.

Procedure:



Results:



Performance:

Statistical Analysis of the QuickVue TLI Lactoferrin Test

Value	QuickVue TLI Lactoferrin vs Lactoferrin Chek®	Active UC* vs IBS and Healthy Persons	Active CD** vs IBS and Healthy Persons
Sensitivity	100.0%	100.0%	100.0%
Specificity	97.5%	100.0%	100.0%
Predictive Positive Value	96.7%	100.0%	100.0%
Predictive Negative Value	100.0%	100.0%	100.0%
Correlation	98.5%	100.0%	100.0%

*Ulcerative colitis

**Crohn's disease

The performance of a rapid Lactoferrin test is dependent on sample collection and handling and the adherence to the Package Insert. Refer to our website at quidel.com for additional performance claims.

1 Sunanda VK, MD, MS, PH et al. Fecal Lactoferrin Is a Sensitive and Specific Marker in Identifying Intestinal Inflammation. Am. J. Gastroenterology: 98 (6) 1309-1314, 2003.