



QuickVue[®] TLI

Technical Bulletin

QuickVue TLI H. pylori Stool Antigen and Campylobacter Universal Reagents

The following reagents are considered to be universal reagents (having the same catalog number) and are common to the QuickVue H. pylori Stool Antigen and Campylobacter products. These reagents are kit-specific but can be interchanged between kit lots. However, we do not recommend the pooling (combining) of different lots of the same reagent due to the inability to track expiration dates of multiple components.

Reagent	Catalog No.	Applicable Kits
Membrane Wash	MW2	20343/20344
Membrane Substrate	MS2	20343/20344

If you have any questions regarding the use of this product, please call Quidel's Technical Support Number 800.874.1517 (in the U.S.) or 858.552.1100, Monday through Friday, from 7:00 a.m. to 5:00 p.m., Pacific Time. If outside the U.S., contact your local distributor or technicalsupport@quidel.com.

You may also visit our website at quidel.com for information on Quidel's line of Rapid Diagnostics, Molecular Diagnostics, Cell Culture and Specialty Products (Bone Health and Autoimmune & Complement). Other product information available on our website includes: CPT codes, CLSI procedure guides, SDS, and Package Inserts.