



Sofia®

Technical Bulletin

Sofia SARS Antigen FIA – Detection of SARS-CoV-2 and SARS-CoV

Sofia SARS Antigen targets a nucleocapsid protein which is present and nearly identical in both SARS-CoV-2, the causative agent of COVID-19, and SARS-CoV, the causative agent of SARS. While the test does not distinguish between the two viruses, SARS-CoV is not actively circulating, according to the CDC, and there have been no known cases of SARS reported anywhere in the world since 2004. This information is available on the CDC website at [cdc.gov/sars](https://www.cdc.gov/sars).

For questions regarding the use of this product, please contact Quidel Technical Support at 1.800.874.1517 (in the U.S.) or technicalsupport@quidel.com. If outside the U.S., further information can be obtained from your distributor, or directly from Quidel at one of the numbers listed below. Reference [quidel.com](https://www.quidel.com) to see more options for Support.

Country	Phone	E-Mail Address
Europe, Middle East and Africa	+353 (91) 412 474 (main) 0 1800 200441 (toll free)	emeatechnicalsupport@quidel.com
Austria	+43 316 231239	
France	0 (805) 371674	
Germany	+49 (0) 7154 1593912	
Netherlands	0 800 0224198	
Switzerland	0 800 554864	
United Kingdom	0 800 3688248	
Italy	+39 (800) 620 549	
North America, Asia-Pacific, Latin America	858.552.1100	technicalsupport@quidel.com
Canada	437.266.1704 (main) 888.415.8764 (toll free)	technicalsupport@quidel.com
China	0400 920 9366 or +86 021 3217 8300	chinatechnicalservice@quidel.com

You may also visit our website at [quidel.com](https://www.quidel.com) for information on Quidel's line of Rapid Diagnostics, Molecular Diagnostics, Cell Culture and Specialty Products (Bone Health and Autoimmune & Complement). Other product information available on our website includes: CPT codes, CLSI procedure guides, SDS, and Package Inserts.