



QuickVue®

Technical Bulletin

Updated Quality Control Recommendations

Effective September 15, 2008, the following new language for quality control frequency recommendations is:

External controls may also be used to demonstrate that the reagents and assay procedure perform properly.

Quidel recommends that positive and negative controls be run once for each untrained operator, once for each new shipment of kits – provided that each different lot received in the shipment is tested – and as deemed additionally necessary by your internal quality control procedures, and in accordance with local, state, and federal regulations or accreditation requirements.

This new language applies to:

- QuickVue Influenza
- QuickVue Influenza A+B
- QuickVue RSV
- QuickVue In-Line® Strep A
- QuickVue Dipstick Strep A

Please contact Quidel Technical Support at 800.874.1517 (in the U.S.), 858.552.1100 (outside the U.S.) or technicalsupport@quidel.com if you have any questions regarding any Quidel product. Our hours of operation are Monday through Friday, 7:00 a.m. to 5:00 p.m. Pacific Time.

You may also visit our website at quidel.com for information on Quidel's line of Rapid Diagnostics, Molecular Diagnostics, Cell Culture and Specialty Diagnostics (Bone Health and Autoimmune & Complement). Other product information available on our website includes: CPT codes, CLSI procedure guides, MSDS, and Package Inserts.