



Technical Bulletin

QuickVue® iFOB: Visibly Bloody Samples

The QuickVue iFOB test is intended for the detection of occult, or "hidden", blood in stool. If the patient has visibly bloody stool, do not use this test.

Please contact Quidel Technical Support at 800-874-1517 (USA only), 858-552-1100, or technicalsupport@quidel.com if you have any questions regarding the QuickVue iFOB test or any other Quidel test. Our hours of operation are Monday-Friday, 7:00 am-5:00 pm Pacific Time.

You may also visit our website at www.quidel.com for information on Quidel's line of Rapid Diagnostics, Bone Health and Autoimmune & Complement product lines. Other product information available on our website includes: CPT codes, CLSI Procedure Guidelines, MSDS and package inserts.