



QuickVue®

Technical Bulletin

QuickVue Influenza and RSV Test Kits: Synagis® Testing

Synagis has been evaluated at a concentration of 1.0 mg/mL and was shown to not interfere with test performance of the following kits:

- QuickVue Influenza (00317)
- QuickVue Influenza A+B (20183)
- QuickVue RSV (20193)
- QuickVue RSV 10 (20222)

Please contact Quidel Technical Support at 800.874.1517 (in the U.S.), 858.552.1100 (outside the U.S.) or technicalsupport@quidel.com if you have any questions regarding the QuickVue Influenza and/or RSV Test, or any Quidel product. Our hours of operation are Monday through Friday, 7:00 a.m. to 5:00 p.m. Pacific Time.

You may also visit our website at quidel.com for information on Quidel's line of Rapid Diagnostics, Molecular Diagnostics, Cell Culture and Specialty Diagnostics (Bone Health and Autoimmune & Complement). Other product information available on our website includes: CPT codes, CLSI procedure guides, MSDS, and Package Inserts.